# Ribospin<sup>™</sup> vRD II

VIRAL RNA PURIFICATION HANDBOOK



### **Customer & Technical Support**

Should you have any further questions, do not hesitate to contact us. We appreciate your comments and advice.

### **Contact Information**

www.geneall.com Tel : 82-2-407-0096 Fax : 82-2-407-0779 E-mail(Order/Sales) : sales@geneall.com E-mail(Tech. Info.) : tech@geneall.com

# Visit GeneAll<sup>®</sup> Community

www.geneall.com www.geneall.co.kr

# Trademarks

Hybrid-Q<sup>™</sup>, Exprep<sup>™</sup>, Exfection<sup>™</sup>, Expin<sup>™</sup>, Exgene<sup>™</sup>, GenEx<sup>™</sup>, DirEx<sup>™</sup>, Hybrid-R<sup>™</sup>, RiboEx<sup>™</sup>, Ribospin<sup>™</sup>, Riboclear<sup>™</sup>, Allspin<sup>™</sup>, RiboSaver<sup>™</sup>, EzClear<sup>™</sup>, EzSep<sup>™</sup>, EzPure<sup>™</sup>, EzPass<sup>™</sup>, AmpONE<sup>™</sup>, AmpMaster<sup>™</sup>, RealAmp<sup>™</sup>, HyperScript<sup>™</sup>, ProtinEx<sup>™</sup>, PAGESTA<sup>™</sup>, STEADi<sup>™</sup>, GENTi<sup>™</sup> are trademarks of GeneAll<sup>®</sup> Biotechnology Co., Ltd.

© 2018 GeneAll Biotechnology, all rights reserved.

This protocol handbook is included in :

GeneAll<sup>®</sup> Ribospin<sup>™</sup> vRD II (322-150, 322-103)

Visit www.geneall.com or www.geneall.co.kr for FAQ, Q&A and more information.

# INDEX

Index	03
Kit Contents	04
Product Specifications	
Quality Control	05
Storage Conditions	
Safety Information	
Preventing RNase Contamination	
Product Description	06
Before Experiment	07
Carrier RNA	
Protocol	08
Troubleshooting Guide	10
Ordering Information	12

# **Kit Contents**

Cat. No.	322-150	322-103	Storage		
Components	Qua	ntity	Storage		
Buffer NVL	l 6 ml	100 ml			
Buffer RB1 (concentrate) *	5 ml	25 ml			
Buffer RBW (concentrate) *	I3 ml	77 ml			
Buffer RNW (concentrate) * <sup>†</sup>	6 ml	34 ml			
Nuclease-free water	l 5 ml	20 ml	Room temperature		
Carrier RNA **	370 µg	2.3 mg	(15~25°C)		
Micro column type S (with collection tube)	50	300			
1.5 ml microcentrifuge tube	50	300			
Protocol Handbook	I.	I.			

\* Before first use, add absolute ethanol (ACS grade or better) into buffer RB I, RBW, and RNW as indicated on the bottle

<sup>+</sup> Contains sodium azide as a preservative

\*\* Refer to page 7 for instruction of Carrier RNA

# **Product Specifications**

Ribospin™ vRD II	
Туре	Spin
Maximum amount of starting samples	Ι 00 <i>μ</i> Ι/prep
Preparation time	$\sim$ I 5 min
Maximum loading volume	750 <i>µ</i> I
Minimum elution volume	20 <i>µ</i> I

### Quality Control

### Storage Conditions

All components in GeneAll<sup>®</sup> Ribospin<sup>™</sup> vRD II are manufactured in strictly clean conditions, and its degree of cleanness is monitored periodically. Quality control is carried out thoroughly from lot to lot, and only the qualified kits are approved to be delivered.

All components of GeneAll<sup>®</sup> Ribospin<sup>TM</sup> vRD II should be stored at room temperature ( $15 \sim 25^{\circ}$ C). After reconstitution of Carrier RNA with Nuclease-free water, the Carrier RNA solution should be stored at -20°C in aliquots for conservation of activity or used immediately for experiments.

During shipment or storage under cool ambient condition, a precipitate can form in buffer NVL. In such a case, heat the bottle to 56°C to dissolve completely. GeneAll<sup>®</sup> Ribospin<sup>TM</sup> vRD II is guaranteed until the expiration data printed on the product box.

Buffer NVL, RBI, and RBW contain an irritant which is harmful when in contact with skin or eyes, or when inhaled or swallowed. Care should be taken when handling such material. Always wear gloves and eye protection, and follow standard safety precautions.

Buffer NVL contains chaotropes, which can form highly reactive compounds when combined with bleach. Do NOT add bleach or acidic solutions directly to the samplepreparation waste.

# Preventing of RNase Contamination

RNase can be introduced accidentally during RNA purification. Wear disposable gloves always, because skin often contains bacteria and molds that can be a source of RNase contamination. Use sterile, disposable plastic wares and automatic pipettes to prevent cross-contamination of RNase from shared equipment.

# Safety Information

### **Product Description**

The GeneAll<sup>®</sup> Ribospin<sup>™</sup> vRD II provides a convenient method for the isolation of RNA and DNA from cell-free fluid, cell-culture supernatant, plasma, serum, swab, urine, and virus-infected samples.

The GeneAll<sup>®</sup> Ribospin<sup>™</sup> vRD II utilizes the glass fiber membrane technology to purify nucleic acid as a sufficient level for downstream analysis instead of conventional alcohol precipitation or phenol/chloroform extraction.

The buffer system of Ribospin<sup>™</sup> vRD II provides the effective binding condition of RNA and DNA to glass fiber membrane and the impurities on the membrane are washed away by two different wash buffers. At last, pure RNA and DNA are eluted in Nuclease-free water. Whole procedure takes only 15 minutes at room temperature and the purified nucleic acid is suitable for PCR, RT-PCR, or any downstream application without further manipulation.

The purified nucleic acid should be treated with care because RNA is very sensitive to contaminants, such as RNases, often found on general labware and dust. To ensure RNA-stability after extraction, it is recommended to store at 4°C for immediate analysis or to freeze at -70°C for long-term storage.

### Before Experiment

Freshly harvested samples should be used or stored immediately for the best result. Starting material, such as plasma or serum, should be stored at  $-70^{\circ}$ C in aliquots for long-term storage.

Repeated freezing and thawing the samples leads to protein precipitation that may cause decreased yield of the extracted viral nucleic acid due to reduction of viral titers in the sample. Besides, the protein precipitant will cause clogging of spin column.

The GeneAll<sup>®</sup> Ribospin<sup>TM</sup> vRD II is designed for extraction of total nucleic acids from samples including virus and host cell. The use of cell-free body fluids is recommended for isolation of viral nucleic acid, and the extraction efficiency can vary depending on the type of virus and sample media.

### Carrier RNA

This kit provides Carrier RNA, which can add at lysis step if required. Provided Carrier RNA can help to improve the binding capacity of mini spin column when viral nucleic acids included in sample are low-copy and protect target nucleic acids from the chance of degradation due to residual RNase activity.

For purification of nucleic acid from very few target molecules in sample, we recommend adding Carrier RNA at lysis step. To obtain a solution of 1  $\mu$ g/ $\mu$ l, add 370  $\mu$ l (Cat. No. 322-150) or 2.3 ml (Cat. No. 322-103) of Nuclease-free water to the tube containing lyophilized Carrier RNA. Dissolve the Carrier RNA thoroughly, divide it into conveniently sized aliquots, and store at -20°C. Do not freeze-thaw the aliquots of Carrier RNA more than 3 times. For one preparation, 7  $\mu$ l of dissolved Carrier RNA is required.

# PROTOCOL FOR **Ribospin<sup>™</sup> vRD II**

#### Before experiment

- Before first use, add absolute ethanol (ACS grade or better) into buffer RB1, RBW and RNW as indicated on the bottle.
- If a precipitate is formed in buffer NVL, heat to  $56^\circ C$  to dissolve before use.
- Prepare an aliquot of Carrier RNA for use on ice (Refer to page 7 for instruction of Carrier RNA).
- ]. Add 300  $\mu$ l of buffer NVL and 7  $\mu$ l of Carrier RNA solution into a 1.5 ml microcentrifuge tube.
- 2. Transfer up to 100  $\mu$ l of sample into the 1.5 ml microcentrifuge tube.

If the sample volume is less than 100  $\mu$ l, adjust the volume to 100  $\mu$ l with 1X PBS. In case of large sample volume, increase the amount of buffer NVL and Carrier RNA solution proportionally.

### 3. Mix thoroughly by vortexing for 10 sec.

For proper lysis, the complete mix of sample and buffer NVL is essential.

#### 4. Incubate the mixture for 10 min at room temperature.

# 5. Add 350 $\mu$ I of buffer RBI to the mixture and mix thoroughly by vortexing for 10 sec.

The volume of buffer RBI can be adjusted in proportion to the volume of lysate.

Normally, the ratio of buffer RB1 to the mixture is 1:1.

Do not centrifuge at this step because nucleic acids can be precipitated through centrifugation.

# 6. Transfer up to 750 $\mu$ l of the mixture to a spin column (Micro column type S, white).

### 7. Centrifuge at $\geq 10,000 \times g$ for 30 sec at room temperature.

Discard the pass-through and reinsert the spin column back into the same tube. If the mixture volume exceeds 750  $\mu$ l, repeat step 6~7 with the remainder.

- 8. Add 500  $\mu$ l of buffer RBW to the spin column.
- 9. Centrifuge at  $\geq 10,000 \times g$  for 30 sec at room temperature.

Discard the pass-through and reinsert the spin column back into the same tube.

- **[0.** Add 500  $\mu$ I of buffer RNW to the spin column.
- | Centrifuge at  $\geq$  10,000 x g for 30 sec at room temperature.

Discard the pass-through and reinsert the spin column back into the same tube.

**12.** Centrifuge at full speed for an additional I min at room temperature to remove residual wash buffer.

#### Transfer the spin column to a new 1.5 ml microcentrifuge tube (provided).

Care must be taken at this step for eliminating the carry-over of buffer RNW that can interfere with downstream reactions.

If a carry-over of buffer RNW still occurs, centrifuge again for 1 min at full speed with the collection tube before transferring to the new 1.5 ml microcentrifuge tube.

# | 3. Add 20~50 $\mu l$ of Nuclease-free water to the center of the membrane in the spin column.

Let it stand for 1 min.

#### 14. Centrifuge at $\geq$ 10,000 x g for 1 min at room temperature.

Purified nucleic acids can be stored at  $4^{\circ}$ C for immediate analysis or can be stored at  $-70^{\circ}$ C for long-term storage.

### Troubleshooting Guide

Facts	Possible Causes	Suggestions
Low yield	Poor quality of starting material	Too old or improperly stored sample often yield degraded DNA. Use fresh sample, if possible. Repeated freezing and thawing the sample should be avoided.
	Low concentration of viral particle in the starting sample	Use more the starting sample. If the amount of sample is more than 300 $\mu$ l, concentrate the volume to 300 $\mu$ l using a microconcentrator.
	Inefficient or insufficient lysis	Be sure to incubate for 10 min at room tem- perature after adding buffer NVL. For proper lysis, the complete mix of sample and buffer NVL is essential.
	Improper elution	Add Nuclease-free water to the center of the spin column membrane and perform incuba- tion for 1 min before centrifugation.
	Precipitation of buffer NVL	Storage at cool ambient temperature may cause precipitation in buffer NVL. For a good result, any precipitate in the buffer should be dissolved by heating the buffer at 56°C or above until it disappears.
	Degradation of RNA	RNase can be introduced during purification of nucleic acid. Be certain not to introduce any RNases during the procedure or later handling. Keep tubes closed whenever pos- sible during the extraction and use RNase- free products with sterile and disposable plastic ware.

### Troubleshooting Guide

Facts	Possible Causes	Suggestions
	Incorrect use of Carrier RNA solution	Add Carrier RNA solution at lysis step. Omis- sion of Carrier RNA may lead to low purifica- tion efficiency.
	Degradation of Carrier RNA	Carrier RNA should be stored at -20°C in ali- quots after reconstitution. Do not freeze-thaw the aliquots of Carrier RNA more than 3 times.
Purified nucleic acid does not perform well	Buffer RBI, RBW, or RNW was prepared incorrectly	Check that the concentrated buffer RBI, RBW, and RNW were diluted with the cor- rect volume of absolute ethanol.
in downstream application	Residual ethanol from buffer RNW remains in eluate	Care must be taken for eliminating the carry- over of buffer RNW before elution step. The membrane of mini spin column should be kept completely dry via additional centrifuga- tion or air-drying.

### Ordering Information

Products	Scale	Size	Cat. No.	Туре	Products	Scale	Size	Cat. No.	Туре
GeneAll® <b>Hybrid</b>	<b>I-Q™</b> fo.	r rapid p	reparation of	plasmid DNA	GeneAll® Exgen	e <sup>™</sup> for is	olation o	f total DNA	
Plasmid Rapidprep		50	100-150				100	105-101	spin /
	mini	200	100-102	spin		mini	250	105-152	vacuum
						MER	26	105-226	spin /
GeneAll® Expre	<b>5<sup>TM</sup></b> for pi	reparatio	on of plasmid	DNA	BIOOD 2V	I*IIdi	100	105-201	vacuum
		50	101-150			MAXI	10	105-310	spin /
	mini	200	101-102	spin /		MAXI	26	105-326	vacuum
		1,000	101-111	vacuum			100	106-101	spin /
Plasmid SV ·		26	101-226		C 11 C /	mini	250	106-152	vacuum
	Midi	50	101-250	spin /	Cell SV		10	106-310	spin /
		100	101-201	vacuum		MAXI	26	106-326	vacuum
							100	108-101	spin /
GeneAll <sup>®</sup> Exfect	ion™					mini	250	108-152	vacuum
for prepa	aration of	transfect	tion-grade pla	ismid DNA			26	108-226	spin /
	mini	50	- 50	spin /	Clinic SV	Midi	100	108-201	vacuum
Plasmid LE		200	- 02	vacuum			10	108-310	spin /
(Low Endotoxin)	Midi	26	-226	spin /		MAXI	26	108-326	vacuum
	THU	100	-20	vacuum	Genomic DNA micr	°0	50	8-050	spin
Plasmid EF	Midi	20	121-220	- spin			100	7- 0	spin /
(Endotoxin Free)	1 IIGI	100	2 -20			mini	250	7- 52	vacuum
							26	117-226	snin /
GeneAll <sup>®</sup> Expin <sup>TI</sup>	<b>Μ</b> for pur	ification	of fragment D	DNA	Plant SV	Midi	100	117-201	vacuum
		50	102-150	spin /			10	117-310	spin /
Gel SV	mini	200	102-102	vacuum		MAXI	26	117-326	vacuum
		50	103-150	spin /	Soil DNA mini	mini	50	114-150	spin
PCR SV	mini	200	103-102	vacuum	Stool DNA mini	mini	50	115-150	spin
		50	113-150	spin /	Viral DNA / RNA	mini	50	128-150	spin
CleanUp SV	mini	200	113-102	vacuum			50	138-150	
		50	2- 50	spin /	FFPE Tissue DNA	mini	250	138-152	spin
Combo GP	mini	200	112-102	vacuum			250	150 152	
					GeneAll® GenEx	тм <sub>for iso</sub>	lation of	total DNA wi	thout spin d
GeneAll <sup>®</sup> Exgene	e <sup>™</sup> for is	olation o	f total DNA				100	220 101	,
		100	04- 0	spin /	GenEx <sup>TM</sup> Blood	Sx	500	220-101	solution
	mini	250	104-152	vacuum	SCHEX BIOOD		100	220 103	solution
		26	104-226	spin /		LA	100	220-301	50101011
Tissue SV	Midi	100	104-201	vacuum	GenEv <sup>TM</sup> Cell	Sx	500	221-101	solution
		10	104-310	spin /	GENEA CEII		100	221-103	colution
	MAXI	26	104-326	vacuum		LX	100	221-301	solutiON
		100	109-101	spin /	ConEv <sup>TM</sup> Tissue	Sx	E00	222-101	solution
	mini	250	109-152	vacuum	Genex Hissue		500	222-105	1 4 <sup>1</sup> .
		200	107-132	acaann		LX	100	222-301	solution

26

MAXI  $\frac{10}{26}$ 

Midi  $\frac{20}{100}$ 

Tissue plus! SV

109-226 spin /

|09-20| vacuum

109-326 vacuum

spin /

109-310

Products	Scale	Size	Cat. No.	Туре

GeneAll<sup>®</sup> GenEx<sup>TM</sup> for isolation of total DNA

	Sx	100	227-101	
GenEx™ Plant	Mx	100	227-201	solution
	Lx	100	227-301	
	Sx	100	228-101	
GenEx <sup>™</sup> Plant plus!	Mx	50	228-250	solution
-	Lx	20	228-320	

GeneAll<sup>®</sup> DirEx<sup>TM</sup> series for preperation of PCR-template without extraction

DirEx™	100	250-101	solution
DirEx <sup>™</sup> <i>Fast-</i> Tissue	96 T	260-011	solution
DirEx <sup>™</sup> <i>Fast</i> -Cultured cell	96 T	260-021	solution
DirEx <sup>™</sup> <i>Fast</i> -Whole blood	96 T	260-03 I	solution
DirEx <sup>™</sup> <i>Fast</i> -Blood stain	96 T	260-041	solution
DirEx <sup>™</sup> <i>Fast</i> -Hair	96 T	260-051	solution
DirEx <sup>™</sup> <i>Fast</i> -Buccal swab	96 T	260-061	solution
DirEx <sup>™</sup> <i>Fast</i> -Cigarette	96 T	260-071	solution

GeneAll®	RNA	series	for preperation	of total	RNA
----------	-----	--------	-----------------	----------	-----

Dibe Ev <sup>TM</sup>	mini	100	301-001	
NIDOEX	(THE)	200	301-002	solution
Hybrid-R <sup>™</sup>	mini	100	305-101	spin
Hybrid-R <sup>™</sup> Blood RNA	٩mini	50	3 5- 50	spin
Hybrid-R <sup>™</sup> miRNA	mini	50	325-150	spin
		100	302-001	a a lu atí a sa
RIDOEX LS	(THE)	200	302-002	solution
Riboclear™	mini	50	303-150	spin
Riboclear <sup>™</sup> plus!	mini	50	3 3- 50	spin
Ribospin™	mini	50	304-150	spin
Dibeenin <sup>™</sup> II	mini	50	3 4- 50	anin
Nibospin II	(THE)	300	3 4- 03	spin
Ribospin <sup>™</sup> vRD	mini	50	302-150	spin
Ribospin <sup>™</sup> vRD <i>plus!</i>	mini	50	3 2- 50	spin
Ribospin <sup>™</sup> vRD II	mini	50	322-150	spin
Ribospin <sup>™</sup> Plant	mini	50	307-150	spin
Ribospin <sup>™</sup> Seed / Fruit	mini	50	317-150	spin
Allspin <sup>™</sup>	mini	50	306-150	spin
RiboSaver™	mini	100	351-001	solution

Products	Scale	Size	Cat. No	<b>)</b> .	Туре
GeneAll® AmpO	NE <sup>™</sup> for	r PCR an	nplificatior	n	
		250 U	501-02	5	
Taq DNA polymeras	se	500 U	501-05	0	(2.5 U/µℓ)
		I,000 U	501-10	0	
		250 U	502-02	5	
lpha-Taq DNA polyme	erase	500 U	502-05	0	(2.5 U/µℓ)
		1,000 U	502-10	0	
-		250 U	504-02	5	
lpha -Pfu DNA polyme	erase	500 U	504-05	0	(2.5 U/µℓ)
		1,000 U	504-10	0	
		250 U	505-02	5	
Fast-Pfu DNA polymerase		500 U	505-05	0	(2.5 U/µℓ)
		1,000 U	505-10	0	
		250 U	531-02	5	
Hotstart Taq DNA		500 U	531-05	0	(2.5 U/µℓ)
polymerase		1,000 U	531-10	0	
		20 µl	521-20	0	h and a literated
T D .	04.1	50 µl	521-50	0	iyopnilized
laq Premix	96 tubes	20 µl	526-20	0	1.2
		50 µl	526-50	0	solution
		20 µl	522-20	0	
	04.1	50 µl	522-50	0	lyophilized
$\alpha$ - Taq Premix	96 tubes	20 μ <b>l</b>	527-20	0	1.2
		50 µl	527-50	0	solution
		20 µl	525-20	0	
HS-Taq Premix	96 tubes	s 50 μl	525-50	0	solution
		20 µl	520-20	0	lyophilized
lpha -Pfu Premix	96 tubes	s 50 μ <b>l</b>	523-50	0	solution
Taq Premix (w/o dye)	96 tubes	s 20 µl	524-20	0	lyophilized
dNTPs mix		500 µl	509-02	0	2.5 mM each
dNTPs set (set of dATP, dCTP, dGTP an	id dTTP)	l ml x 4 tubes	509-04	10	100 mM

Products	Scale	Size	Cat. No.	Туре
----------	-------	------	----------	------

### GeneAll<sup>®</sup> AmpMaster<sup>™</sup> for PCR amplification

Taq Master mix	0.5 ml x 2 tubes	541-010	solution
	0.5 ml x 10 tubes	541-050	solution
lpha -Taq Master mix	0.5 ml x 2 tubes	542-010	solution
	0.5 ml x 10 tubes	542-050	solution
HS-Taq Master mix	0.5 ml x 2 tubes	545-010	solution
	0.5 ml x 10 tubes	545-050	solution
lpha-Pfu Master mix	0.5 ml x 2 tubes	543-010	solution
	0.5 ml x 10 tubes	543-050	solution

### GeneAll<sup>®</sup> HyperScript<sup>™</sup> for Reverse Transcription

Reverse Transcript	ase 10,000 U	601-100	solution
RT Master mix	0.5  ml  imes 2  tubes	601-710	solution
RT Master mix with oligo (dT) <sub>20</sub>	$0.5 \ \mathrm{ml}  imes 2 \ \mathrm{tubes}$	601-730	solution
RT Master mix with random hexamer	$0.5 \ \mathrm{ml}  imes 2 \ \mathrm{tubes}$	601-740	solution
RT Premix	96 tubes, 20 µl	601-602	solution
RT Premix with oligo (dT) <sub>20</sub>	96 tubes, 20 μ <b>l</b>	601-632	solution
RT Premix with random hexamer	96 tubes, 20 µl	601-642	solution
One-step RT-PCR Master mix	$0.5 \ \mathrm{ml}  imes 2 \ \mathrm{tubes}$	602-110	solution
One-step RT-PCR Premix	96 tubes, 20 μ <b>l</b>	602-102	solution
First strand Synthesis Kit	50 reaction	605-005	solution
ZymAll <sup>™</sup> RNase Inhibitor	10,000 U	605-010	solution
ZymAll <sup>™</sup> RNase Inhibitor	4,000 U	605-004	solution

### GeneAll<sup>®</sup> RealAmp<sup>™</sup> for qPCR amplification

SYBR qPCR Master	200 rxn	20 <i>µl</i>	801-020	
mix (2X, Low ROX)	500 rxn	20 <i>µ</i> l	801-050	solution
SYBR qPCR Master	200 rxn	20 <i>µl</i>	801-021	achtion
mix (2X, High ROX)	500 rxn	20 µl	801-051	SOlution

Products	Size	Cat. No.
----------	------	----------

### GeneAll<sup>®</sup> Protein series

ProtinEx <sup>™</sup> Animal cell / tissue	100 ml	701-001	solution
PAGESTA <sup>™</sup> Reducing 5X SDS-PAGE I mI > Sample Buffer	< 10 tubes	751-001	solution

### $\mathsf{GeneAll}^{\circledast}\operatorname{\mathsf{STEAD}}\!\!i^{ imes}$ for automatic nucleic acid puritication

12 Instrument		GST012	system
24 Instrument		GST024	system
Genomic DNA Cell / Tissue	96	401-104	kit
Genomic DNA Blood	96	402-105	kit
Bacteria DNA	96	403-106	kit
Total RNA	96	404-304	kit
Viral DNA / RNA	96	405-322	kit
CFC Seed DNA / RNA	96	406-C02	kit
Genomic DNA Plant	96	407-107	kit
Soil DNA	96	407-108	kit

### NOTE



www.geneall.com

GeneAll Bldg., 303-7 Dongnam-ro, Songpa-gu, Seoul, South Korea 05729 E-mail : sales@geneall.com

Tel : 82-2-407-0096 Fax : 82-2-407-0779

©2018 GeneAll Biotechnology, All rights reserved