

Cat. No. 307-150

Ribospin™ Plant

PLANT TOTAL RNA PURIFICATION HANDBOOK

Customer & Technical Support

Do not hesitate to ask us any question.

We thank you for any comment or advice.

Contact us at

www.geneall.com

Tel : 82-2-407-0096

Fax : 82-2-407-0779

E-mail(Order/Sales) : sales@geneall.com

E-mail(Tech. Info.) : tech@geneall.com

Visit GeneAll® Community

www.geneall.com

www.geneall.co.kr

Trademarks

Hybrid-Q™, Exprep™, Exfection™, Expin™, Exgene™, GenEx™, DirEx™, Hybrid-R™, RiboEx™, Ribospin™, Ribodear™, Allspin™, RiboSaver™, EzClear™, EzSep™, EzPure™, EzPass™, AmpONE™, AmpMaster™, RealAmp™, HyperScript™, ProtinEx™, PAGESTA™, STEADi™ are trademarks of GeneAll Biotechnology co., ltd.

© 2016 GeneAll Biotechnology, all right reserved.

This protocol handbook is included in :

GeneAll® Ribospin™ Plant (307-150)

Visit www.geneall.com or www.geneall.co.kr for FAQ, QnA and more information.

Homogenization

Grind ~100 mg/prep plant samples in liquid nitrogen.

Transfer the powder into a 1.5 ml microcentrifuge tube.

Lysis step

Add 350 μ l of buffer RPL.

Incubate the mixture for 3 min at RT.

EzPure™ filter step

Transfer the lysate to a EzPure™ filter and centrifuge at $\geq 10,000 \times g$ for 30 sec.

Transfer the supernatant into a 1.5 ml microcentrifuge tube.

RNA binding step

Add 1 volume of 70% ethanol to the supernatant and mix well.

Apply the mixture into a mini spin column and centrifuge at $\geq 10,000 \times g$ for 30 sec.

DNase I treatment step

Add 500 μ l of buffer RBW to the mini spin column and centrifuge at $\geq 10,000 \times g$ for 30 sec.

Apply the DNase I mixture into a mini spin column.

Incubate the mixture for 10 min at RT.

Washing step

Add 500 μ l of buffer RBW to the mini spin column and incubate for 2 min and centrifuge at $\geq 10,000 \times g$ for 30 sec.


Add 500 μ l of buffer RNW to the mini spin column and centrifuge at $\geq 10,000 \times g$ for 30 sec (twice).

Centrifuge at $\geq 10,000 \times g$ for an additional 1 min.

RNA elution

Add ~50 μ l of Nuclease-free water to the center of the membrane.

Centrifuge at $\geq 10,000 \times g$ for 1 min.



INDEX

	Brief Protocol	03
	Index	05
	Kit Contents	06
	Materials Not Provided	
	Quality Control	07
	Storage Conditions	
	Precautions	
	Preventing RNase Contamination	
	Product Disclaimer	
	Product Specifications	08
	Typical Yields	
	Product Description	09
	Protocol for Ribospin™ Plant	10
	Troubleshooting Guide	13
	Appendix 1	15
	Appendix 2	16
	Appendix 3	17
	Ordering Information	19

KIT CONTENTS

Ribospin™ Plant kit

Components	Quantity	Storage
Buffer RPL	25 ml	Room temperature
Buffer REL	25 ml	
Buffer RBW (concentrate) *	27 ml	
Buffer RNW (concentrate) * †	12 ml	
Nuclease-free water	15 ml	
Buffer DRB	5 ml	
GeneAll® EzPure™ filter (yellow) (with collection tube)	50	
GeneAll® Column type W (blue ring) (with collection tube)	50	
1.5 ml microcentrifuge tube	100	
DNase I	120 ul	

* Before using for the first time, add absolute ethanol (ACS grade or better) into buffer RBW and RNW as indicated on the bottle.

† Contains sodium azide as a preservative

Materials Not Provided

Reagent

- 70% ethanol, ACS grade or better
- Absolute ethanol, ACS grade or better

Disposable material

- RNase-free pipet tips
- Disposable gloves

Equipment

- Equipment for disrupting plant tissue
- Microcentrifuge

Quality Control

Ribospin™ Plant kit is manufactured in strictly clean condition, and its degree of cleanness is monitored periodically. For consistency of product, the quality certification process is carried out from lot to lot thoroughly and only the qualified is approved to be delivered.

Storage Conditions

Ribospin™ Plant kit, except DNase I, should be stored at room temperature (15~25°C). DNase I should be stored at -20°C. All components are stable for 1 year.

Precautions

The buffers included in Ribospin™ Plant kit contain irritant which is harmful when in contact with skin or eyes, or when inhaled or swallowed. Care should be taken during handling. Always wear gloves and eye protector, and follow standard safety precautions. In case of contact, wash immediately with plenty of water and seek medical advice. Buffer RPL, REL, and RBW contain chaotropes. It can form highly reactive compounds when combined with bleach. Do NOT add bleach or acidic solutions directly to the sample-preparation waste.

Preventing RNase Contamination

RNase can be introduced accidentally into a RNA preparation. Wear disposable gloves always, because skin often contains bacteria that can be a source of RNase. Use sterile, disposable plasticwares and automatic pipettes reserved for RNA work to prevent cross-contamination with RNase on shared equipment.

Product Disclaimer

Ribospin™ Plant kit is for research use only. This kit is not to be used for any other clinical test such as diagnostic, prognostic, therapeutic, etc.

Product Specifications

Specification	Ribospin™ Plant
Type	Spin
Maximum amount of starting samples	~ 100 mg plant tissue
Maximum loading volume of spin column	~700 ul
Minimum elution volume	30 ul
Maximum binding capacity	~ 100 ug

Typical Yields

	Sample type	Amount of starting material	Typical yield
Leaf	<i>Pinus densiflora</i> (Pine)	100 mg	2.7 ug
	<i>Cucumis sativus</i> L. (Cucumber)	100 mg	50 ug
	<i>Zea mays</i> (Corn)	100 mg	11 ug
	<i>Capsicum annuum</i> (Red pepper)	100 mg	22 ug
	<i>Lycopersicon esculentum</i> (Tomato)	50 mg	13 ug
	<i>Lactuca sativa</i> (Lettuce)	100 mg	29 ug
	<i>Citrus grandis</i> Osbek (Satsuma)	100 mg	4.6 ug
	<i>Diospyros kaki</i> (Persimmon)	100 mg	16 ug
	<i>Crassula ovata</i> (Crassula)	100 mg	3 ug
	<i>Nicotiana tabacum</i> (Tobacco)	50 mg	13 ug
Root	<i>Allium cepa</i> (Onion)	100 mg	8 ug
	<i>Plantago asiatica</i> (Plantain)	50 mg	2.5 ug
	<i>Nicotiana tabacum</i> (Tobacco)	50 mg	5.3 ug
Fruit	<i>Citrus grandis</i> Osbek (Satsuma)	50 mg	1.1 ug
Germ bud	<i>Allium cepa</i> (Onion)	100 mg	9 ug

Product Description

Ribospin™ Plant kit is specially designed for purification of total RNA from various plant tissues such as leaves, stems, roots and picky plant samples. This kit provides the optimized buffer and spin column, which is effective at removing polysaccharides and polyphenolic compounds and isolating intact plant RNA. All components of Ribospin™ Plant are ready for use, so any further preparation for experiment is not required.

The procedure of Ribospin™ Plant kit begins with the disruption of sample in liquid nitrogen using mortar and pestle. The disrupted sample can be lysed in buffer RPL or REL. In most case, buffer RPL is the best buffer for lysis. However in some plant samples, solidification of lysate can be occurred with buffer RPL due to endosperm of seed or peculiar metabolites, and this can be avoided by using buffer REL as alternative for buffer RPL.

Most impurities except RNA in the lysate are eliminated by filtration through EzPure™ filter, and then the passed-through lysate is mixed with ethanol to adjust binding condition. Total RNA including a little impurity is bound to the membrane of spin column type W while the mixture is passing through. Survived genomic DNA can be exterminated by on-column DNase I treatment at this step. After a series of washing step using buffer RBW and RNW, plant total RNA is eluted by Nuclease-free water.

Whole procedure of Ribospin™ Plant takes only 25 minutes. The purified RNA is suitable for cDNA synthesis, RT-PCR, Northern blotting, and other analytical procedure.

PROTOCOL FOR Ribospin™ Plant

Before experiment

- Before using for the first time, add absolute ethanol (ACS grade or better) into buffer RBW and RNW as indicated on the bottle.
- Prepare DNase I reaction mixture just before step 12.
 - √ Prepare aliquot DNase I and thaw on ice.
 - √ Mix 2 ul DNase I with 70 ul Buffer DRB.

- 1. Prepare plant tissue sample up to 100 mg, then grind the sample to a fine powder using a mortar and pestle with liquid nitrogen and transfer the grinded sample into a 1.5 ml microcentrifuge tube (not provided).**
- 2. Add 350 ul of buffer RPL to the 1.5 ml microcentrifuge tube and vortex vigorously.**

In case of solidification of the lysate in buffer RPL, use buffer REL instead of buffer RPL.
- 3. Incubate 3 min at room temperature.**
- 4. Transfer the lysate to a EzPure™ filter.**

Through this step, large cell debris and most of genomic DNAs are filtered on the EzPure™ filter and small pellet of cell debris will be formed at the bottom of the collection tube.
- 5. Centrifuge at $\geq 10,000 \times g$ for 30 seconds at room temperature.**
- 6. Transfer the supernatant to a new 1.5 ml microcentrifuge tube (provided).**

Be careful not to disturb the pellet at the bottom of the collection tube.
- 7. Add 1 volume (usually 350 ul) of 70% EtOH to the tube containing supernatant, and mix well by pipetting or inverting.**

Do not centrifuge at this step.

8. Apply the mixture to a mini spin column (type W, blue ring).

9. Centrifuge at $\geq 10,000 \times g$ for 30 seconds at room temperature.

Discard the pass-through and reinsert the mini spin column back into the same tube.

10. Add 500 ul of buffer RBW to the mini spin column.

11. Centrifuge at $\geq 10,000 \times g$ for 30 seconds at room temperature.

Discard the pass-through and reinsert the mini spin column back into the same tube.

12. Apply 70 ul of DNase I reaction mixture to the center of the mini spin column.

Incubate at the room temperature for 10 minutes.

To make DNase I reaction mixture, mix 2 ul DNase I with 70 ul Buffer DRB. DNase I is sensitive to physical damage and thus do not mix vigorously. If you want to DNase I treatment in RNA eluate, skip step 12 and 13 and refer to Appendix I "DNase I treatment in eluate".

13. Add 500 ul of buffer RBW to the mini spin column and stand for 2 minutes.

Buffer RBW inactivates DNase I and wash out the components of DNase I reaction buffer.

14. Centrifuge at $\geq 10,000 \times g$ for 30 seconds at room temperature.

Discard the pass-through and reinsert the mini spin column back into the same tube.

15. Add 500 ul of buffer RNW to the mini spin column.

16. Centrifuge at $\geq 10,000 \times g$ for 30 seconds at room temperature.

Discard the passed-through and reinsert the mini spin column back into the same tube.

17. Repeat step 15~16.

18. Centrifuge at $\geq 10,000 \times g$ for an additional 1 minute at room temperature to remove residual wash buffer. Transfer the mini spin column to a new 1.5 ml microcentrifuge tube (provided).

Residual ethanol may interfere with downstream reactions. Care must be taken at this step for eliminating the carryover of buffer RNW.

19. Add 50 ul of Nuclease-free water to the center of the membrane in the mini spin column.

To increase the RNA concentration, reduce the elution volume to 30 ul.

20. Centrifuge at $\geq 10,000 \times g$ for 1 minute at room temperature.

Purified RNA can be stored at 4°C for immediate analysis and can be stored at -70°C for long term storage.

The purified RNA is free of DNA and proteins, and A_{260}/A_{280} will be between 1.8 and 2.2.

Troubleshooting Guide

Facts	Possible Causes	Suggestions
Low yield of RNA	Sample not disrupted completely.	Insufficient disruption can lead to decrease the yield of total RNA. Confirm the completely disrupted sample in liquid nitrogen and transfer the disrupted sample in a 1.5 ml microcentrifuge tube.
	Too much starting sample	Overloading can decrease the yield of total RNA. Reduce the amount of starting sample.
	Poor quality of starting material	Process the sample immediately after harvest. To process later, freeze the sample rapidly in liquid nitrogen.
	Too low RNA mass in samples	Especially, some plant samples have low RNA content. To increase the RNA concentration, reduce the elution volume up to 30 ul or increase the amount of starting sample up to 100 mg per prep.
RNA degradation	Sample manipulated too much before process	Process the sample immediately after harvest. To process later, freeze the sample rapidly in liquid nitrogen.
	Improper storage of RNA	Store isolated RNA at -70°C, Do not store at -20°C.
	Reagent or disposable is not RNase-free	Make sure to use RNase-free products only.
EzPure™ filter clogging	Sample not disrupted completely	Insufficient disruption can clog the EzPure™ filter and to decrease the yield of total RNA. Confirm the complete disruption of the sample in liquid nitrogen.

Troubleshooting Guide

Facts	Possible Causes	Suggestions
DNA contamination of RNA eluate	Too much starting sample	Too much starting sample may leave lots of DNA fragments on the membrane over the activity of DNase I. Reduce starting sample used.
	Sample has high DNA mass	Some plant samples have high DNA content. In this case, some DNA can be eluted at RNA elution step. Reducing the amount of sample can reduce the genomic DNA contamination or refer to the appendix I 'DNase I treatment in eluate'.
	DNase I not active	For prolonged activity, aliquot the DNase I into small portion. Do not freezing and thawing the aliquots several times.
	Incorrect DNase I reaction treatment	Add DNase I reaction mixture to the center of the mini spin column membrane.
Eluate does not perform well in downstream application	Residual ethanol remains in eluate	To remove any residual ethanol included in buffer RNW from mini spin column membrane, centrifuge again for complete removal of ethanol.
	Buffer RBW and RNW used in wrong order	Ensure that buffer RBW and RNW are used in correct order. If used in the wrong order, wash the spin column with buffer RNW finally.

APPENDIX I. DNase I treatment in eluate

Appendix I describe how to use the DNase I (included in this kit) to eliminate contaminating DNA in RNA eluate. For samples containing high DNA contents, this method is strongly recommended. This procedure is more efficient than on-column DNase I treatment.

Protocol

1. Prepare the mixture as below in a 1.5 ml microcentrifuge tube.

50 ul RNA eluate

5 ul Buffer DRB

1 ul DNase I

2. Incubate the mixture for 10 minutes at room temperature.

3. Re-elution of RNA.

Follow 3-1 or 3-2

DNase I treated RNA can be applied to RNA clean up kit (Riboclear™ Cat. No. 303-150).

We strongly recommend using Riboclear™ kit for RNA clean up. Because ethanol precipitation and heat inactivation, usually used for DNase I inactivation, can damage the RNA.

3-1 Follow Riboclear™ protocol

3-2 Heat inactivation

1. Add 1 ul of 0.5 M EDTA per 100 ul eluate.

2. Heat inactivate at 75°C for 10 minutes.

APPENDIX 2.

Confirmation of RNA yield and purity by UV absorbance

Concentration of RNA

The concentration of RNA can be determined by the absorbance at 260 nm using spectrophotometer. For the convenient measurement, we recommend using the NanoDrop® which can reduce your RNA sample and time. If not, you need to dilute the RNA samples to measure the concentration through traditional spectrophotometer. The value of A_{260} should be between 0.15 and 1.00. Be sure to calibrate the spectrophotometer with the same solution used for dilution. An absorbance of 1 at 260 nm is correspond to about 40 ug RNA/ml at a neutral pH. Therefore, the concentration of RNA was calculated by the formula shown below.

$$A_{260} \times \text{dilution factor} \times 40 = \text{RNA ug/ml}$$

Purity of RNA

To confirm the RNA purity, you should read the ratio of A_{260}/A_{280} . Pure RNA is in the range of 1.8~2.2.

APPENDIX 3. Formaldehyde agarose gel electrophoresis (Denaturing gel method)

A denaturing agarose gel is routinely used for the assessment of the quality of an RNA preparation. After preparation, RNA forms secondary structure via intramolecular base pairing. Therefore, it is very difficult to get the exact result of electrophoresis because of migrating inaccuracy. However, the denaturing gel denatures the secondary structure of RNA and makes an accurate migration.

To confirm the RNA band, the gel should be transferred to a UV transilluminator after electrophoresis. Mainly, two RNA bands are shown. In case of animal sample, the 28S and 18S rRNA bands are confirmed on the gel. If they are intact, the RNA bands should be sharp and the intensity of upper band should be about twice that of the lower band.

Prepare the denaturing gel

1. Put 1g agarose in 72 ml water and heat to dissolve thoroughly.
2. Cool to 60°C.
3. Add 10 ml of 10X MOPS buffer, 18 ml of 37% formaldehyde, and 1 ul of a 10 mg/ml ethidium bromide (EtBr).
4. Mix well then pour the gel into the gel tray and cool to solidify it.
5. Transfer the solidified gel from tray to tank, and add enough 1X MOPS running buffer to cover the gel.

Prepare the RNA sample

1. Make the mixture.
 - x ul RNA (up to 20 ug)
 - 2 ul 10X MOPS electrophoresis buffer
 - 4 ul Formaldehyde
 - 10 ul Formamide
2. Incubate the mixture for 15 minutes at 65°C.
3. Chill the sample for 5 minutes on ice.
4. Add 2 ul of 10X formaldehyde gel-loading dye to the mixture.
5. Load the mixture in a denaturing gel which is covered with a sufficient 1X MOPS electrophoresis buffer.
6. Run the gel and confirm the RNA band on transilluminator.
Occasionally, gel destaining may be needed to increase the visibility of the bands of RNA in dH₂O for several hours.

Composition of buffers

- 10X MOPS buffer

0.2 M MOPS

20 mM Sodium acetate

10 mM EDTA

pH to 7.0 with NaOH

- 10X formaldehyde gel-loading dye

50% Glycerol

10 mM EDTA

0.25% (w/v) Bromophenol blue

0.25% (w/v) Xylene cyanol FF

* **Caution**

When working with these chemicals, always use gloves and eye protector to avoid contact with skin and cloth. Especially, formaldehyde and ethidium bromide (EtBr) should be handled in a fume hood.

Ordering Information

Products	Scale	Size	Cat. No.	Type
----------	-------	------	----------	------

GeneAll® Hybrid-Q™ for rapid preparation of plasmid DNA

Plasmid Rapidprep	mini	50	100-150	spin
		200	100-102	

GeneAll® Exprep™ for preparation of plasmid DNA

	mini	50	101-150	spin /	
		200	101-102	vacuum	
Plasmid SV	Midi	26	101-226	spin /	
		50	101-250		vacuum
		100	101-201		

GeneAll® Exfection™ for preparation of transfection-grade plasmid DNA

Plasmid LE (Low Endotoxin)	mini	50	111-150	spin /
		200	111-102	vacuum
	Midi	26	111-226	spin /
		100	111-201	vacuum
Plasmid EF (Endotoxin Free)	Midi	20	121-220	spin
		100	121-201	

GeneAll® Expin™ for purification of fragment DNA

Gel SV	mini	50	102-150	spin /
		200	102-102	vacuum
PCR SV	mini	50	103-150	spin /
		200	103-102	vacuum
CleanUp SV	mini	50	113-150	spin /
		200	113-102	vacuum
Combo GP	mini	50	112-150	spin /
		200	112-102	vacuum

GeneAll® Exgene™ for isolation of total DNA

Tissue SV	mini	100	104-101	spin /
		250	104-152	vacuum
	Midi	26	104-226	spin /
		100	104-201	vacuum
	MAXI	10	104-310	spin /
		26	104-326	vacuum
Tissue plus! SV	mini	100	109-101	spin /
		250	109-152	vacuum
	Midi	26	109-226	spin /
		100	109-201	vacuum
	MAXI	10	109-310	spin /
		26	109-326	vacuum

Products	Scale	Size	Cat. No.	Type
----------	-------	------	----------	------

GeneAll® Exgene™ for isolation of total DNA

Blood SV	mini	100	105-101	spin /
		250	105-152	vacuum
	Midi	26	105-226	spin /
100		105-201	vacuum	
	MAXI	10	105-310	spin /
		26	105-326	vacuum
Cell SV	mini	100	106-101	spin /
		250	106-152	vacuum
	MAXI	10	106-310	spin /
		26	106-326	vacuum
Clinic SV	mini	100	108-101	spin /
		250	108-152	vacuum
	Midi	26	108-226	spin /
		100	108-201	vacuum
	MAXI	10	108-310	spin /
		26	108-326	vacuum
Genomic DNA micro		50	118-050	spin
Plant SV	mini	100	117-101	spin /
		250	117-152	vacuum
	Midi	26	117-226	spin /
		100	117-201	vacuum
	MAXI	10	117-310	spin /
		26	117-326	vacuum
Soil DNA mini	mini	50	114-150	spin
Stool DNA mini	mini	50	115-150	spin
Viral DNA / RNA	mini	50	128-150	spin
FFPE Tissue DNA	mini	50	138-150	spin
		250	138-152	

GeneAll® GenEx™ for isolation of total DNA without spin column

GenEx™ Blood	Sx	100	220-101	solution
		500	220-105	
	Lx	100	220-301	solution
GenEx™ Cell	Sx	100	221-101	solution
		500	221-105	
	Lx	100	221-301	solution
GenEx™ Tissue	Sx	100	222-101	solution
		500	222-105	
	Lx	100	222-301	solution

Products	Scale	Size	Cat. No.	Type
GeneAll® GenEx™ for isolation of total DNA				
GenEx™ Plant	Sx	100	227-101	solution
	Mx	100	227-201	
	Lx	100	227-301	
GenEx™ Plant plus!	Sx	100	228-101	solution
	Mx	50	228-250	
	Lx	20	228-320	

GeneAll® DirEx™ series
for preparation of PCR-template without extraction

Products	Scale	Size	Cat. No.	Type
DirEx™		100	250-101	solution
DirEx™ Fast-Tissue	96 T		260-011	solution
DirEx™ Fast-Cultured cell	96 T		260-021	solution
DirEx™ Fast-Whole blood	96 T		260-031	solution
DirEx™ Fast-Blood stain	96 T		260-041	solution
DirEx™ Fast-Hair	96 T		260-051	solution
DirEx™ Fast-Buccal swab	96 T		260-061	solution
DirEx™ Fast-Cigarette	96 T		260-071	solution

GeneAll® RNA series for preparation of total RNA

Products	Scale	Size	Cat. No.	Type
RiboEx™	mini	100	301-001	solution
		200	301-002	
Hybrid-R™	mini	100	305-101	spin
Hybrid-R™ Blood RNA mini	50		315-150	spin
Hybrid-R™ miRNA	mini	50	325-150	spin
RiboEx™ LS	mini	100	302-001	solution
		200	302-002	
Riboclear™	mini	50	303-150	spin
Riboclear™ plus!	mini	50	313-150	spin
Ribospin™	mini	50	304-150	spin
Ribospin™ II	mini	50	314-150	spin
		300	314-103	
Ribospin™ vRD	mini	50	302-150	spin
Ribospin™ vRD plus!	mini	50	312-150	spin
Ribospin™ vRD II	mini	50	322-150	spin
Ribospin™ Plant	mini	50	307-150	spin
Ribospin™ Seed / Fruit	mini	50	317-150	spin
Allspin™	mini	50	306-150	spin
RiboSaver™	mini	100	351-001	solution

Products	Scale	Size	Cat. No.	Type
GeneAll® AmpONE™ for PCR amplification				
Taq DNA polymerase		250 U	501-025	(2.5 U/μl)
		500 U	501-050	
		1,000 U	501-100	
α-Taq DNA polymerase		250 U	502-025	(2.5 U/μl)
		500 U	502-050	
		1,000 U	502-100	
α-Pfu DNA polymerase		250 U	504-025	(2.5 U/μl)
		500 U	504-050	
		1,000 U	504-100	
Fast-Pfu DNA polymerase		250 U	505-025	(2.5 U/μl)
		500 U	505-050	
		1,000 U	505-100	
Hotstart Taq DNA polymerase		250 U	531-025	(2.5 U/μl)
		500 U	531-050	
		1,000 U	531-100	
Taq Premix	96 tubes	20 μl	521-200	lyophilized
		50 μl	521-500	
		20 μl	526-200	
		50 μl	526-500	
α-Taq Premix	96 tubes	20 μl	522-200	lyophilized
		50 μl	522-500	
		20 μl	527-200	
		50 μl	527-500	
HS-Taq Premix	96 tubes	20 μl	525-200	solution
		50 μl	525-500	
		20 μl	520-200	
α-Pfu Premix	96 tubes	50 μl	523-500	solution
Taq Premix (w/o dye)	96 tubes	20 μl	524-200	lyophilized
dNTPs mix		500 μl	509-020	2.5 mM each
dNTPs set (set of dATP, dCTP, dGTP and dTTP)		1 ml x 4 tubes	509-040	100 mM

Products	Scale	Size	Cat. No.	Type
----------	-------	------	----------	------

GeneAll® AmpMaster™ for PCR amplification

Taq Master mix	0.5 ml × 2 tubes	541-010	solution
	0.5 ml × 10 tubes	541-050	solution
α-Taq Master mix	0.5 ml × 2 tubes	542-010	solution
	0.5 ml × 10 tubes	542-050	solution
HS-Taq Master mix	0.5 ml × 2 tubes	545-010	solution
	0.5 ml × 10 tubes	545-050	solution
α-Pfu Master mix	0.5 ml × 2 tubes	543-010	solution
	0.5 ml × 10 tubes	543-050	solution

GeneAll® HyperScript™ for Reverse Transcription

Reverse Transcriptase	10,000 U	601-100	solution
RT Master mix	0.5 ml × 2 tubes	601-710	solution
RT Master mix with oligo (dT) ₂₀	0.5 ml × 2 tubes	601-730	solution
RT Master mix with random hexamer	0.5 ml × 2 tubes	601-740	solution
RT Premix	96 tubes, 20 μl	601-602	solution
RT Premix with oligo (dT) ₂₀	96 tubes, 20 μl	601-632	solution
RT Premix with random hexamer	96 tubes, 20 μl	601-642	solution
One-step RT-PCR Master mix	0.5 ml × 2 tubes	602-110	solution
One-step RT-PCR Premix	96 tubes, 20 μl	602-102	solution
First strand Synthesis Kit	50 reaction	605-005	solution
ZymAll™ RNase Inhibitor	10,000 U	605-010	solution
ZymAll™ RNase Inhibitor	4,000 U	605-004	solution

GeneAll® RealAmp™ for qPCR amplification

SYBR qPCR Master mix (2X, Low ROX)	200 rxn 20 μl	801-020	solution
	500 rxn 20 μl	801-050	
SYBR qPCR Master mix (2X, High ROX)	200 rxn 20 μl	801-021	solution
	500 rxn 20 μl	801-051	

Products	Size	Cat. No.
----------	------	----------

GeneAll® Protein series

ProteinEx™ Animal cell / tissue	100 ml	701-001	solution
PAGESTA™ Reducing 5X SDS-PAGE Sample Buffer	1 ml × 10 tubes	751-001	solution

GeneAll® STEADi™ for automatic nucleic acid purification

STEADi™ 12 Instrument		GST012
STEADi™ 24 Instrument		GST024
STEADi™ Genomic DNA Cell / Tissue Kit	96	401-104
STEADi™ Genomic DNA Blood Kit	96	402-105
STEADi™ Bacteria DNA Kit	96	403-106
STEADi™ Total RNA Kit	96	404-304
STEADi™ Viral DNA / RNA Kit	96	405-322
STEADi™ CFC Seed DNA / RNA Kit	96	406-C02

Note

Note



www.geneall.com

GeneAll Bldg., 303-7 Dongnam-ro,
Songpa-gu, Seoul, South Korea | 38-859
E-mail : sales@geneall.com

Tel : 82-2-407-0096
Fax : 82-2-407-0779

©2016 GeneAll Biotechnology, All right reserved

Edited by **BnP**
Designed by **KDY**
2016.08