## Amplicon Genotyping and **Methylation** Analysis by HRM

Jonathan Nery **Applications Scientist** 

© 2010 Illumina, Inc. All rights reserved. Illumina, illuminaDx, Solexa, Making Sense Out of Life, Oligator, Sentrix, GoldenGate, GoldenGate Indexing, DASL, BeadArray, Array of Arrays, Infinium, BeadXpress, VeraCode, IntelliHyb, iSelect, CSPro, GenomeStudio, Genetic Energy, HiSeg, and HiScan are registered trademarks or trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.



#### Introduction

- Requirements for High Resolution Melt
  - Chemistry
  - Instrumentation
  - Software
- Applications
  - Genotyping
    - Class IV SNPs
    - Open platform
  - DNA Methylation Analysis



#### **Requirements for HRM**

#### Chemistry

- Saturating dsDNA binding dyes
- Non-inhibitory to PCR
- SYTO<sup>®</sup> 9 (Life Tech.), LCGreen<sup>™</sup> (Idaho Tech.), EvaGreen<sup>™</sup> (Biotium Inc.)

#### Instrumentation

- Precise temperature control
- Minimum of well-to-well thermal or optical non uniformity
- Fast data acquisition rate for more data points

#### Software

- Specialized normalization algorithms
- Difference plots



## **Double Stranded DNA Binding Dyes**

Non-Saturating dsDNA Binding Dye - SYBR<sup>®</sup> Green I



illumina

Saturating dsDNA Binding Dye

- SYTO<sup>®</sup> 9

#### **HRM Compatible Real-Time PCR Instruments**



Life Technologies 7500Fast/7900Fast/Viia 7



Illumina Eco™



Qiagen Rotor-Gene Q



Bio-Rad CFX96/CFX384



Roche LightCycler® 480

Images not to scale



#### **Eco Thermal System**



- Temperature uniformity and accuracy of < ±0.1°C</p>
- Ramp rate of > 5°C/sec
- Temperature resolution of 0.1°C
- No temperature shift or calibration required for HRM



#### Eco Software – Melt Curve Normalization

- Raw fluorescence pre- and post- melt normalized to 100% and 0%
- Aligns curves for better differentiation





#### **Eco Software – Difference Plot View**

- Alternate view of melt curve to maximize differences
- Subtract a reference curve from each curve





## **Applications – Amplicon Genotyping**

- Three possible genotypes
  - Homozygous Wild-type
  - Homozygous Mutant
  - Heterozygous



Each will produce a unique melting profile that can be differentiated by HRM



#### **PCR of Homozygous Samples**





#### **Melting Profile of Homozygous Samples**



Homozygous Samples Differentiated by Shift in Melting Temperature



11

#### **PCR of Heterozygous Sample**



#### **Melting Profile of Heterozygous Samples**



Heterozygous Sample Differentiated by Curve Shape



13

## Amplicon $T_m$ Shifts for the 4 SNP Classes

SNP Class	Base Change	Frequency in Human Genome	Approximate T <sub>m</sub> Melt Curve Shift
1	C/T and G/A	0.67	> 0.5°C
2	C/A and G/T	0.18	0.4 to 0.5°C / 0.2 to 0.5°C
3	C/G	0.09	0.2 to 0.4°C
4	A/T	0.07	< 0.2°C

- Amplicon T<sub>m</sub> shifts vary from large to very small depending on identity of polymorphism
- Class 4 SNPs are the rarest and most difficult SNP type to discriminate by T<sub>m</sub>
- HRM requires highly optimized assays and precise instrumentation to differentiate sequence changes



#### **Single Nucleotide Resolution**

#### Primer and Amplicon Sequences (60 and 61 bp):

5'-	CACCTCACGCAGCACTTACCAA	
5'-	CACCTCACGCAGCACTTACCAACTACTCAT <b>a</b> CAGACTCATTCACCTCACCATGTCACTCGC	76.64°C
5 <b>'</b> -	CACCTCACGCAGCACTTACCAACTACTCAT <b>t</b> CAGACTCATTCACCTCACCATGTCACTCGC	76.82°C
5'-	CACCTCACGCAGCACTTACCAACTACTCAT <b>C</b> CAGACTCATTCACCTCACCATGTCACTCGC	77.47°C
5'-	CACCTCACGCAGCACTTACCAACTACTCATgCAGACTCATTCACCTCACC	77.70°C
5'-	CACCTCACGCAGCACTTACCAACTACTCAT : CAGACTCATTCACCTCACCATGTCACTCGC	76.95°C
	<b>AAGTGGAGTGGTACAGTGAGCG</b> $-5'$	



#### **Normalized Melt and Difference Plot**



Five clear melt profiles with good separation



A/T polymorphism within KIAA1683 (rs8110972)

Primer and Amplicon Sequences (48 bp):

5' - GGTGACAGCCATGTCTACA	
5'- GGTGACAGCCATGTCTACAGGCACA <b>a</b> ACACGTGAGGI	GGCTTGTCCCC 81.17°C
5'- GGTGACAGCCATGTCTACAGGCACAtaCACGTGAGGI	GGCTTGTCCCC 80.76°C
ACTCCA	CCGAACAGGGG -5 '

▶ 8 human genomic DNA samples – 3 technical replicates of each



Normalized Melt Normalized Fluorescence minus Reference Normalized Fluorescence -5 -10 **BLUE** = TT PINK= AA **GREEN** = AT -15 0 -Temperature (°C) Temperature (°C)

**Difference Plot** 

#### illumina

 A/T polymorphism within ATP-binding cassette, sub-family C (CFTR/MRP), member 12 (ABCC12) gene (rs6500305)

Primer and Amplicon Sequences (75 bp):

5′-	CCAGGCCCTGCATGGA	
5′-	CCAGGCCCTGCATGGAGGAGGTGATGTGGG <b>t</b> GAACCAGGGTGACCGGCTGACATTCTCCACCTTCTTGAGCTCCT	84.08°C
5′-	CCAGGCCCTGCATGGAGGAGGTGATGTGGG <mark>a</mark> GAACCAGGGTGACCGGCTGACATTCTCCACCTTCTTGAGCTCCT	83.95°C
	TAAGAGGTGGAAGAACTCGAGGA	-5′

45 human genomic DNA samples



Normalized Melt Normalized Fluorescence minus Reference Normalized Fluorescence -5 -10 **BLUE** = TT PINK= AA -15 **GREEN** = AT Temperature (°C) Temperature (°C)

**Difference Plot** 

#### illumina

#### **Available HRM Master Mixes**

Four commercially available HRM mixes evaluated on Eco system

- BioRad SsoFast EvaGreen Supermix
- Life Technologies EXPRESS SYBR<sup>®</sup> GreenER<sup>™</sup>
- Life Technologies MeltDoctor™ HRM Master Mix
- Roche LightCycler <sup>®</sup> 480 HRM
- Genotyping of a Class II SNP
  - A/C polymorphism in solute carrier family 28 (sodium-coupled nucleoside transporter), member 2 (SLC28A2) gene (rs1060896)

#### Primer and Amplicon Sequences (65 bp):

5'- AAAGCAAGAAGTTTCTGCAAAACA	
5'- AAAGCAAGAAGTTTCTGCAAAACACACGCCAG <mark>a</mark> TTGTTCAAGAAGATCCTGTTGGGCCTGTTGTG	77.12°C
5'- AAAGCAAGAAGTTTCTGCAAAACACACGCCAG <mark>c</mark> TTGTTCAAGAAGATCCTGTTGGGCCTGTTGTG	77.93°
TAGGACAACCCGGACAACAC -5'	

45 human genomic DNA samples



#### **Normalized Melt Curves**



Three genotypes clearly resolved by all mixes



#### **Difference Plots**



Wild Type (AA) sample used a reference for difference plot



23

### **Application – DNA Methylation Analysis**

 DNA methylation of (CpGs) is a key epigenetic mechanism regulating gene expression





## **Methylation Analysis by HRM**



Samples differentiated by differences in melting temperature



#### **HRM can be Sensitive and Semi-Quantitative**

- Bisulfite converted methylated and unmethylated DNA (Qiagen) mixed in different ratios to form a 'standard curve'
- Primers designed according to recommendations from Wojdacz et al (2009)
  - Limited number of CpGs included in primers to overcome PCR bias
  - 60°C annealing temperature

#### **HRM can be Sensitive and Semi-Quantitative**





#### **HRM can be Sensitive and Semi-Quantitative**



Both assays able to detect 1% methylated DNA in a background of unmethylated DNA



# Methylation Analysis of Tumor, Pseudotumor and Normal Tissue

**Tumor sample**: expected to be less methylated than the normal adjacent tissue **Pseudotumor sample**: expected to have same methylation level as the normal adjacent tissue





#### Summary

- Requires high performance instrumentation
- Enables cost effective high value analysis
  - Genotyping
  - Methylation







illumına<sup>®</sup>

30

#### Acknowledgements

The Eco team



McGill University

- Moshe Szyf, Ph.D.
- Barbara Stefanska, Ph.D
- Aurélie Bouzelmat

