Ver 4.4





#### **Customer & Technical Support**

Should you have any further questions, do not hesitate to contact us.

We appreciate your comments and advice.

#### **Contact Information**

www.geneall.com

Tel: 82-2-407-0096

Fax: 82-2-407-0779

E-mail(Order/Sales): sales@geneall.com

E-mail(Tech. Info.): tech@geneall.com

#### Visit GeneAll® Community

www.geneall.com

www.geneall.co.kr

#### **Trademarks**

Hybrid- $Q^{TM}$ , Exprep<sup>TM</sup>, Exfection<sup>TM</sup>, Expin<sup>TM</sup>, Exgene<sup>TM</sup>, GenEx<sup>TM</sup>, DirEx<sup>TM</sup>, Hybrid-R<sup>TM</sup>, RiboEx<sup>TM</sup>, RiboSpin<sup>TM</sup>, RiboSaver<sup>TM</sup>, EzClear<sup>TM</sup>, EzSep<sup>TM</sup>, EzPure<sup>TM</sup>, EzPass<sup>TM</sup>, AmpONE<sup>TM</sup>, AmpMaster<sup>TM</sup>, RealAmp<sup>TM</sup>, HyperScript<sup>TM</sup>, ProtinEx<sup>TM</sup>, PAGESTA<sup>TM</sup>, STEADi<sup>TM</sup> are trademarks of GeneAll Biotechnology Co., Ltd.

© 2017 GeneAll Biotechnology, all rights reserved.


This protocol handbook is included in :

GeneAll® Expin<sup>™</sup> Gel SV (102-150, 102-102) GeneAll® Expin<sup>™</sup> PCR SV (103-150, 103-102) GeneAll® Expin<sup>™</sup> CleanUp SV (113-150, 113-102) GeneAll® Expin<sup>™</sup> Combo GP (112-150, 112-102)

Visit www.geneall.com or www.geneall.co.kr for FAQ, Q&A and more information.

## **Brief Protocol**

## **Expin<sup>™</sup> Gel SV**



- I. Excise the DNA band of interest
- 2. Add 3 vol. (µI) Buffer GB to I vol. (mg) gel
- 3. Incubate at 50°C until the gel is completely melted



- 4. Apply the mixture into mini column
- 5. Centrifuge for I min



- 6. Add 700 μl Buffer NW
- 7. Apply the mixture into mini column
- 8. Centrifuge for 30 sec




9. Additional centrifuge for 1 min



- I 0. Apply 50 μl Buffer EB
- II. Incubate for I min at RT
- 12. Centrifuge for I min

## **Expin™ PCR SV**



I. Add 5 vol. Buffer PB to I vol. sample



- 2. Apply the mixture into mini column
- 3. Centrifuge for 30 sec



- 4. Add 700 μl Buffer NW
- 5. Apply the mixture into mini column
- 6. Centrifuge for 30 sec




7. Additional centrifuge for I min



- 8. Apply 50  $\mu$ l Buffer EB
- 9. Incubate for I min at RT
- 10. Centrifuge for 1 min

## Expin<sup>™</sup> Clean up SV



I. Add I0 vol. Buffer NR to I vol. sample



- 2. Apply the mixture into mini column
- 3. Centrifuge for 30 sec



- 4. Add 700  $\mu$ l Buffer NW
- 5. Apply the mixture into mini column



7. Additional centrifuge for I min



- 8. Apply 50  $\mu$ l Buffer EB
- 9. Incubate for I min at RT
- 10. Centrifuge for I min

## 

## GeneAll<sup>®</sup> Expin<sup>™</sup> Protocols

Brief Protocol	3
Index	7
Kit Contents	8
Product Specifications	8
Quality Control	9
Storage Conditions	9
Safety Information	9
Product Disclaimer	9
Product Description	10
- Applicable range in length of DNA	- 11
- General Considerations	12
- Binding	12
- Washing	13
- Elution	13
- Correlation between the elution volume and the recovery	14
- Electrophoresis analysis	15
PROTOCOLS FOR	
$Expin^TM Gel  SV  Protocol$	16
$Expin^TM$ PCR SV Protocol	22
$Expin^TM$ CleanUp SV Protocol	26
Troubleshooting Guide	30
Ordering Information	32

## **Kit Contents**

	Ge	I SV	PCR	SV
Cat. No.	102-150	102-102	103-150	103-102
Size	mini	mini	mini	mini
No. of preparation	50	200	50	200
Column Type D (mini) (with collection tube)	50	200	50	200
Buffer GB	60 ml	120 ml x 2	-	-
Buffer PB	-	-	30 ml	120 ml
Buffer NW (concentrate) * †	12 ml	50 ml	12 ml	50 ml
Buffer EB **	15 ml	30 ml	15 ml	30 ml
Protocol Handbook	1	1	1	1

	Clean	Up SV	Comb	o GP
Cat. No.	113-150	113-102	112-150	112-102
Size	mini	mini	mini	mini
No. of preparation	50	200	50	200
Column Type D (mini) (with collection tube)	50	200	50	200
Buffer GB	-	-	60 ml	120 ml x 2
Buffer PB	-	-	30 ml	120 ml
Buffer NR	30 ml	120 ml	-	-
Buffer NW (concentrate) * †	12 ml	50 ml	12 ml	50 ml
Buffer EB **	15 ml	30 ml	15 ml	30 ml
Protocol Handbook	1	1	1	1

<sup>\*</sup> Before first use, add absolute ethanol (ACS grade or better) into Buffer NW as indicated on the bottle

## **Product Specifications**

GeneAll® Expin™	Gel SV	PCR SV	CleanUp SV
Туре	spin/mini	spin/mini	spin/mini
Standard sample size	200 mg of gel	100 $\mu$ l of PCR	$50  \mu \text{I}$ of enzyme
	(400 mg max.)	reactions	mixtures
Recovered DNA size	80 bp~10 kb	100 bp∼10 kb	40 bp~10 kb
Typical yields	≤85%	≤95%	≤95%
Binding capacity	10 μg	10 μg	10 μg
Preparation time	≥15 min	≥6 min	≥6 min
Maximum loading volume of mini column	750 <i>μ</i> l	750 μl	$750~\mu$ l
Minimum elution volume	30 <i>μ</i> l	30 $\mu$ l	30 <i>μ</i> Ι

<sup>&</sup>lt;sup>†</sup> Contains sodium azide as a preservative

<sup>\*\* 10</sup> mM TrisCl, pH 8.5

## Quality Control

All components in GeneAll® Expin<sup>TM</sup> series are manufactured in strictly clean conditions, and its degree of cleanness is monitored periodically. Quality control is carried out thoroughly from lot to lot, and only the qualified kits are approved to be delivered.

## Storage Conditions

All components of GeneAll® Expin<sup>TM</sup> series should be stored at room temperature ( $15\sim25$ °C). It should be protected from exposure to direct sunlight.

During shipment or storage under cool ambient condition, a precipitate can form in Buffer GB, PB or NR. In such a case, heat the bottle to  $50^{\circ}$ C to dissolve completely. Using precipitated buffers will lead to poor DNA recovery. GeneAll<sup>®</sup> Expin<sup>TM</sup> series are guaranteed until the expiration date printed on the product box.

## Safety Information

The buffers included in the GeneAll® Expin<sup>TM</sup> series contain irritants which is harmful when in contact with skin or eyes, or when inhaled or swallowed. Care should be taken when handling such materials. Always wear gloves and eye protection, and follow standard safety precautions.

Buffer GB, PB and NR contain chaotropes, which can form highly reactive compounds when combined with bleach. Do NOT add bleach or acidic solutions directly to the sample-preparation waste.

## Product Disclaimer

GeneAll® Expin $^{TM}$  series is for research use only, not for use in diagnostic procedure.



#### **Product Description**

GeneAll® Expin™ series provide reliable and fast methods for the purification of fragment DNA from agarose gel and PCR/enzyme reaction mixtures. GeneAll® Expin™ series consist of;

#### GeneAll<sup>®</sup> Expin<sup>™</sup> Gel SV (102-150, 102-102)

GeneAll® Expin™ Gel SV kit is designed for a fast and efficient isolation of 80 bp to 10 kb of DNA fragments from standard or low-melting agarose gel in TAE or TBE buffer system. Purified DNA can be directly used in ligation, labelling, sequencing and many other downstream application without further manipulation.

#### **GeneAll® Expin<sup>™</sup> PCR SV (103-150, 103-102)**

GeneAll® Expin<sup>TM</sup> PCR SV kit provides a simple and rapid method to purify PCR products or other enzymatic reactions in just 6 minutes. Up to 10  $\mu$ g of pure DNA (100 bp~10 kb) can be obtained, and this purified DNA can be directly used in cloning, sequencing and many other application. GeneAll® Expin<sup>TM</sup> PCR SV procedures remove the DNA fragment smaller than 100 bp, resulting in removal of primers and primer dimers in PCR products.


#### GeneAll® Expin<sup>™</sup> CleanUp SV (113-150, 113-102)

GeneAll® Expin<sup>TM</sup> CleanUp SV kit provides a simple and rapid method to purify 40 bp~10 kb DNA from enzymatic reactions in just 6 minutes. Up to 10  $\mu$ g of pure DNA which is at least 40 bp but less than 10 kb in length can be obtained using this kit and the purified DNA can be directly used for sequencing, cloning and other routine applications without further manipulation.

#### GeneAll® Expin<sup>™</sup> Combo GP (112-150, 112-102)

The combination set of GeneAll® Expin™ Gel SV and Expin™ PCR SV.

### - Applicable range in length of DNA



#### - General Considerations

GeneAll® Expin™ series takes advantage of silica membrane and spin/vacuum column technology to recover DNA fragments. Under high salt conditions, DNA binds to silica membrane and impurities pass through membrane into a collection tube. The membranes are washed with an ethanol-containing buffer to remove any traces of proteins, salts, remnants, of agarose and other enzymatic reaction components. Finally pure DNA is released into a clean collection tube with water or low ionic strength buffer.

#### - Binding

The basic principle which DNA binds to silica membrane is identical in all GeneAll® Expin<sup>TM</sup> series. Binding buffers (Buffer GB, PB and NR) make the optimal binding condition in each specific applications. Buffer GB in Expin™ Gel SV kit is composed to dissolve standard agarose gel as well as low melt agarose gel in addition to adjusting a binding condition. Usually low melt agarose gel results in a better recovery yield. DNA binds to silica membrane at lower pH than pH 7.5. The components and pH of starting sample can alter the pH of the mixture with binding buffer, especially in extraction of DNA from agarose gel. Buffer GB contains pH indicator in order to check this alteration of binding condition. If the color (yellow) of binding mixture turns to brown or purple after addition of Buffer GB, it means that the pH of binding mixtures is higher than the optimal, and it can be easily adjusted with small volume of sodium acetate before proceeding with the protocol. The indicator dye is completely removed during subsequent washing steps and does not interfere the downstream applications.

#### - Washing


Any unwanted oligos and impurities, such as salts, proteins, nucleotides, agarose, dyes and detergents will not be bound but be passed through the silica membrane. A minute impurities, such as salts, are quantitatively washed away with Buffer NW which contains ethanol. The quality of DNA can be slightly increased with the repeat of washing. Any residual ethanol should be removed completely with an additional centrifugation because the residual ethanol in eluate may interfere some subsequent applications.

#### - Elution

DNA is released under the condition of low salts and neutral or weakly alkaline pH (7.0<pH<9.5). Although Buffer EB (10 mM TrisCl, pH 8.5), TE, or distilled water can be used for elution, it should be considered that EDTA in buffer TE may interfere the subsequent reactions and low pH (<7.0) of distilled water can reduce DNA recovery. Because water does not have any buffering agents the eluate in water should be stored under -20°C not to degrade.

The minimum elution volume is 30  $\mu$ l and lower volume will decrease the yield significantly. It is important for optimal elution to apply the elution buffer to the center of the membrane, because the membrane should be covered completely by eluent for an optimal recovery. Up to 200  $\mu$ l of elution buffer can be applied and it results in low concentration of DNA. Higher concentrated DNA will be obtained with lower elution volume, and maximum yield can be obtained by larger elution volume. The yield with large fragments (>5 kb) can be increased slightly by using pre-warmed (70 °C) elution buffer. Incubation for 1 minute after addition of eluent may increase the efficiency of elution.

#### - Correlation between the elution volume and the recovery




Correlation between the elution volume and the recover rates. Upper panel; GeneAll® Expin™ PCR SV kit. Lower panel; GeneAll® Expin™ Gel SV kit. 3 µg of 494 bp PCR products were purified and eluted with the indicated volumes of Buffer EB. Elution volume lower than 30 µl causes significant loss of DNA. 1/10 volume of eluate purified with Expin™ Gel SV kit was resolved on 1% agarose gel; lower right. (1/20 for Expin™ PCR SV kit, upper right)

	Recovery rates (%)					
DNA size (bp)	Gel SV	PCR SV	CleanUp SV			
60	39	0	63			
120	71	78	80			
200	76	83	84			
800	84	94	94			
1800	82	91	93			
4300	78	85	88			
8700	73	76	79			

Average recovery rates of GeneAll® Expin<sup>TM</sup> series with various sizes of DNA. 3 μg of starting sample was purified and eluted with 50 μl of Buffer EB. Optional steps were not performed and mini columns were incubated for 1 minute after addition of Buffer EB.

#### - Electrophoresis analysis



Quantities of extracted 4.5 kb DNA fragment corresponding to 1/5 of the DNA obtained by purification from 0.5  $\mu$ g starting DNA with a recovery of 90%. Sample were run on 1% TAE agarose gel. M; Lambda-BstP1 marker Lane1; Total amount before extraction (0.5  $\mu$ g) Lane2; 1/5 amount after extraction [90 ng compared to known amount (10~150 ng) DNA] \* Total obtained amount of DNA=90 x 5=450 ng approximately (90%)



PCR products of several sizes (70, 176, 757 and 1487 bp from left to right) were purified using GeneAll® Expin<sup>™</sup> CleanUp SV kit. The band intensity of 1/5 amount of 20 µl PCR products (lane 1, 5, 9, 13) was almost identical with that of 1/5 amount of 50 µl eluate. (Lane 2~4, 6~8, 10~12 and 14~16, tripli-cate) M; 1 kb ladder marker

#### Expin<sup>™</sup> Gel SV Protocol

#### Before experiment

- \* Before first use, add absolute ethanol (ACS grade or better) to Buffer NW as indicated on the bottle.
- \* All centrifugation should be carried out at  $10,000 \times g$  above (>12,000 rpm) at room temperature in a microcentrifuge.
- \* Prepare water bath or heating block to 50°C.
- \* All solutions should be equilibrated at room temperature before procedures.
- \* For large fragments (>5 kb), pre-warm Buffer EB at 70°C.
- \* If a precipitate is formed in Buffer GB, heat at 50 °C to dissolve before use.

#### Spin/Vacuum Methods for Gel Extraction

When using this kit, one of the two methods can be chosen to purify DNA. DNA fragment can be purified using centrifugation to pull the mixture through the mini column. Alternatively, vacuum can be used to force the mixture out through the column.

#### A. Centrifugation Protocol

I. Excise the DNA band of interest using a ethanol-cleaned razor blade or scalpel on a transilluminator.

Minimize gel volume by cutting the gel slice as small as possible.

Use of long wave length transilluminator and short handling time will lead to better quality of DNA. It can be critical in some experiments, such as ligation.

2. Weigh the gel slice in a 1.5 ml microcentrifuge tube. Add 3 volumes  $(\mu I)$  of Buffer GB to I volume (mg) of gel.

For example, add 300  $\mu$ l of Buffer GB to each 100 mg of agarose gel slice. For > 1.5% agarose gel, add 5 volumes of Buffer GB.

3. Incubate at  $50^{\circ}$ C until the agarose gel is completely melted (5~10 min).

To help the efficient dissolving of gel, vortex the tube every  $2\sim3$  min during the incubation.

4. After the slice has dissolved completely, check that the color of the mixture is yellow (similar to Buffer GB).

If the color of the mixture becomes brown or purple, add  $10 \mu l$  of 3 M sodium acetate, pH 5.0, and mix. The color of the mixture will turn to yellow. Adjusting of pH is not needed when the change of color is due to the ingredient of loading dye (eg. bromophenol blue, xylene cyanol).

5. (Optional:) Add I gel volume of isopropanol to the sample and vortex to mix.

For 100 mg of gel volume, add 100  $\mu$ l of isopropanol.

Do NOT centrifuge at this step.

This step is required to increase the recover yields of DNA fragments <200 bp or >5 kb. For the DNA fragments between 200 bp and 5 kb, it has little effect on the recovery.

6. Transfer the mixture to a Column Type D (mini). Centrifuge for I min at  $10,000 \times g$  above (>12,000 rpm). Discard the pass-through and re-insert the mini column into the collection tube.

If the mixture volume is larger than 700  $\mu$ l, apply the mixture twice; apply 700  $\mu$ l of the mixture, spin down, discard the pass-through, re-insert empty collection tube, and repeat the step again until all of the mixture has been applied to the mini column.

7. (Optional:) Apply 500  $\mu$ l of Buffer GB to the mini column. Centrifuge for 30 sec at 10,000 x g above (>12,000 rpm). Discard the pass-through and re-insert the mini column into the collection tube.

This step is for further complete removal of any traces of agarose and required only for direct use of purified DNA for very sensitive applications, such as in vitro transcription. Usually this step can be skipped for automatic sequencing or ligation.

8. Add 700  $\mu$ l of Buffer NW to the mini column. Centrifuge for 30 sec at 10,000 x g above (>12,000 rpm). Discard the pass-through and re-insert the mini column into the collection tube.

If the purified DNA will be used for salt sensitive applications, let the mini column stand for 5 min after addition of Buffer NW, making some amount of wash buffer flow through the column by gravity before centrifugation.

 Centrifuge for an additional I min at full speed to remove residual wash buffer. Transfer the mini column to a new 1.5 ml microcentrifuge tube.

If residual ethanol remains in the mini column, centrifuge again for an additional I min at full speed before transferring to a new 1.5 ml microcentrifuge tube. Residual ethanol from Buffer NW can inhibit subsequent enzymatic reaction.

10. Apply 50  $\mu$ I of Buffer EB or ddH<sub>2</sub>O to the center of the membrane in the mini column, let stand for I min and centrifuge for I min at 10,000 x g above (>12,000 rpm).

Ensure that the Buffer EB or distilled water is dispensed directly onto the center of column membrane for optimal elution of DNA.

To obtain more concentrated DNA solution, apply 30  $\mu$ l of elution buffer, but the volume lower than 30  $\mu$ l will decrease the yield significantly. Up to 200  $\mu$ l of elution buffer can be applied to mini column, and it will reduce the concentration of DNA.

For larger fragment (>5 kb), use pre-warmed (70°C) elution buffer for best efficiency.

For long-term storage, eluting in Buffer EB (10 mM TrisCl, pH 8.5) or TE, pH 8.0 and storing under -20 $^{\circ}$ C is recommended. When using water for elution, check the pH of water (>pH 7.0) before elution.

#### **B.** Vacuum Protocol

The vacuum pressure should be in the range of the list below. Lower vacuum pressure may reduce DNA yield and purity.

Most commercial vacuum manifold with luer connectors can be used with this protocol.

I. Excise the DNA band of interest using a ethanol-cleaned razor blade or scalpel on a transilluminator.

Minimize gel volume by cutting the gel slice as small as possible.

Use of long wave length transilluminator and short handling time will lead to better quality of DNA. It can be critical in some experiments, such as ligation.

2. Weigh the gel slice in a microcentrifuge tube. Add 3 volumes ( $\mu$ I) of Buffer GB to I volume (mg) of gel.

For example, add 300  $\mu$ l of Buffer GB to each 100 mg of agarose gel slice. For >1.5% agarose gel, add 5 volumes of Buffer GB.

3. Incubate at  $50^{\circ}$ C until the agarose gel is completely melted (5~10 min).

To help the efficient dissolving of gel, vortex the tube every  $2\sim3$  min during the incubation.

4. After the slice has dissolved completely, check that the color of the mixture is yellow (similar to Buffer GB).

If the color of the mixture becomes brown or purple, add  $10 \,\mu$ I of 3 M sodium acetate, pH 5.0, and mix. The color of the mixture will turn to yellow. Adjusting of pH is not needed when the change of color is due to the ingredient of loading dye (eg. bromophenol blue, xylene cyanol)

5. (Optional:) Add I gel volume of isopropanol to the sample and vortex to mix.

For 100 mg of gel volume, add 100  $\mu$ l of isopropanol. Do NOT centrifuge at this step.

This step is required to increase the recover yields of DNA fragments <200 bp or >5 kb. For the DNA fragments between 200 bp and 5 kb, it has little effect on the recovery.

6. Attach a Column Type D (mini) to a port of the vacuum manifold tightly.

Most commercial vacuum manifold with luer connectors can be adopted to this protocol.

- 7. Transfer the mixture to the mini column by pipetting. Switch on vacuum source to draw the solution through the mini column. When all liquid has been pulled through the mini column, release the vacuum. Repeat this step until all mixture have applied to the mini column.
- 8. (Optional:) Apply 500  $\mu$ l of Buffer GB to the mini column and switch on vacuum source. When all liquid has been pulled through the mini column, release the vacuum.

This step is for further complete removal of any traces of agarose and required only for direct use of purified DNA for very sensitive applications, such as in vitro transcription. Usually this step can be skipped for automatic sequencing.

9. Apply 700  $\mu$ l of Buffer NW and switch on vacuum source. When all liquid has been pulled through the mini column, release the vacuum. Transfer the mini column to a empty collection tube (provided).

If the purified DNA will be used for salt sensitive applications, let the mini column stand for 5 minutes after addition of Buffer NW before applying vacuum.

[10. Centrifuge for an additional I min at full speed to remove residual wash buffer. Transfer the mini column to a new 1.5 ml microcentrifuge tube.

If the mini column has Buffer NW associated with it, centrifuge again for additional I min at full speed before transferring to a new 1.5 ml microcentrifuge tube. Residual ethanol from Buffer NW can inhibit subsequent enzymatic reaction.

## | | . Apply 50 $\mu$ l of Buffer EB or ddH<sub>2</sub>O to the center of the membrane in the mini column, let stand for 1 min and centrifuge for 1 min at 10,000 x g above (>12,000 rpm).

Ensure that the Buffer EB or distilled water is dispensed directly onto the center of column membrane for optimal elution of DNA.

To obtain more concentrated DNA solution, apply 30  $\mu$ l of elution buffer, but the volume lower than 30  $\mu$ l will decrease the yield significantly. Up to 200  $\mu$ l of elution buffer can be applied to mini column, and it will reduce the concentration of DNA.

For larger fragment (>5 kb), use pre-warmed (70°C) elution buffer for best efficiency.

For long-term storage, eluting in Buffer EB (10 mM TrisCl, pH 8.5) or TE, pH 8.0 and storing under -20°C is recommended. When using water for elution, check the pH of water (>pH 7.0) before elution.

#### **Expin<sup>™</sup> PCR SV Protocol**

#### Before experiment

- \* Before first use, add absolute ethanol (ACS grade or better) to Buffer NW as indicated on the bottle.
- \* All centrifugation should be carried out at  $10,000 \times g$  above (>12,000 rpm) at room temperature in a microcentrifuge.
- \* All solutions should be equilibrated at room temperature before procedures.
- \* For large fragments (>5 kb), pre-warm Buffer EB at 70°C.

#### Spin/Vacuum Methods for PCR Purification

When using this kit, one of the two methods can be chosen to purify DNA. DNA fragment can be purified using centrifugation to pull the mixture through the mini column. Alternatively, vacuum can be used to force the mixture through the column.

#### A. Centrifugation Protocol

I. Add 5 volumes of Buffer PB to I volume of the sample and mix. Transfer the mixture to a Column Type D (mini).

For  $100 \,\mu l$  reaction, add  $500 \,\mu l$  of Buffer PB. It is not necessary to remove mineral oil.

- 2. Centrifuge for 30 sec at  $10,000 \times g$  above (>12,000 rpm). Discard the pass-through and re-insert the mini column back into the same tube.
- 3. Apply 700  $\mu$ l of Buffer NW and centrifuge for 30 sec at 10,000 x g above (>12,000 rpm). Discard the pass-through and re-insert the mini column back into the collection tube.

 Centrifuge for an additional I min at full speed to remove residual wash buffer. Transfer the mini column to a new 1.5 ml microcentrifuge tube.

If the mini column has Buffer NW associated with it, centrifuge again for additional I min at full speed before transferring to the new 1.5 ml microcentrifuge tube. Residual ethanol from Buffer NW can inhibit subsequent enzymatic reaction.

5. Apply 50  $\mu$ I of Buffer EB or ddH<sub>2</sub>O to the center of the membrane in the mini column, let stand for I min and centrifuge for I min at 10,000 x g above (>12,000 rpm).

Ensure that the Buffer EB or distilled water is dispensed directly onto the center of column membrane for optimal elution of DNA.

To obtain more concentrated DNA solution, apply 30  $\mu$ l of elution buffer, but the volume lower than 30  $\mu$ l will decrease the yield significantly. Up to 200  $\mu$ l of elution buffer can be applied to mini column, and it will reduce the concentration of DNA.

For larger fragment (>5 kb), use pre-warmed ( $70^{\circ}$ C) elution buffer for best efficiency.

For long-term storage, eluting in Buffer EB (10 mM TrisCl, pH 8.5) or TE, pH 8.0 and storing under -20°C is recommended. When using water for elution, check the pH of water (>pH 7.0) before elution.

#### **B.** Vacuum Protocol

The vacuum pressure should be in the range of the list below. Lower vacuum pressure may reduce DNA yield and purity.

Most commercial vacuum manifold with luer connectors can be used with this protocol.

- I. Attach a Column Type D (mini) to a port of the vacuum manifold tightly.
- Add 5 volumes of Buffer PB to I volume of the sample and mix. Transfer the mixture to the mini column by pipetting.

For 100  $\mu$ l reaction, add 500  $\mu$ l of Buffer PB.

It is not necessary to remove mineral oil.

- 3. Switch on vacuum source to draw the solution through the mini column. When all liquid has been pulled through the mini column, release the vacuum.
- 4. Apply 700  $\mu$ l of Buffer NW and switch on vacuum source. When all liquid has been pulled through the mini column, release the vacuum. Transfer the mini column to a collection tube (provided).
- 5. Centrifuge for I min at full speed to remove residual wash buffer. Transfer the mini column to a new 1.5 ml microcentrifuge tube.

If the mini column has Buffer NW associated with it, centrifuge again for additional I min at full speed before transferring to the new 1.5 ml microcentrifuge tube. Residual ethanol from Buffer NW can inhibit subsequent enzymatic reaction.

# 6. Apply 50 $\mu$ l of Buffer EB or ddH<sub>2</sub>O to the center of the membrane in the mini column, let stand for 1 min and centrifuge for 1 min at 10,000 x g above (>12,000 rpm).

Ensure that the Buffer EB or distilled water is dispensed directly onto the center of column membrane for optimal elution of DNA.

To obtain more concentrated DNA solution, apply 30  $\mu$ l of elution buffer, but the volume lower than 30  $\mu$ l will decrease the yield significantly. Up to 200  $\mu$ l of elution buffer can be applied to mini column, and it will reduce the concentration of DNA.

For larger fragment (>5 kb), use pre-warmed (70°C) elution buffer for best efficiency.

For long-term storage, eluting in Buffer EB (10 mM TrisCl, pH 8.5) or TE, pH 8.0 and storing under -20°C is recommended. When using water for elution, check the pH of water (>pH 7.0) before elution.

#### Expin<sup>™</sup> CleanUp SV Protocol

#### Before experiment

- \* Before first use, add absolute ethanol (ACS grade or better) to Buffer NW as indicated on the bottle.
- \* All centrifugation should be carried out at  $10,000 \times g$  above (>12,000 rpm) at room temperature in a microcentrifuge.
- \* All solutions should be equilibrated at room temperature before procedures.
- \* For large fragments (>5 kb), pre-warm Buffer EB at 70°C.

#### Spin/Vacuum Methods for DNA Clean-up

When using this kit, one of the two methods can be chosen to purify DNA. DNA fragment can be purified using centrifugation to pull the mixture through the mini column. Alternatively, vacuum can be used to force the mixture through the column.

#### A. Centrifugation Protocol

I. Add 10 volumes of Buffer NR to I volume of the sample and mix. Transfer the mixture to a Column Type D (mini).

For 50  $\mu$ l reaction, add 500  $\mu$ l of Buffer NR. If the length of DNA is longer than 100 bp, add 5 volumes of Buffer NR.

- 2. Centrifuge for 30 sec at 10,000 x g above (>12,000 rpm). Discard the pass-through and re-insert the mini column back into the collection tube.
- 3. Apply 700  $\mu$ l of Buffer NW and centrifuge for 30 sec at 10,000 x g above (>12,000 rpm). Discard the pass-through and re-insert the mini column back into the collection tube.

4. Centrifuge for an additional I min at full speed to remove residual wash buffer. Transfer the mini column to a new 1.5 ml microcentrifuge tube.

If the mini column has Buffer NW associated with it, centrifuge again for additional I min at full speed before transferring to the new 1.5 ml microcentrifuge tube. Residual ethanol from Buffer NW can inhibit subsequent enzymatic reaction.

5. Apply 50  $\mu$ I of Buffer EB or ddH<sub>2</sub>O to the center of the membrane in the mini column, let stand for I min and centrifuge for I min at 10,000 x g above (>12,000 rpm).

Ensure that the Buffer EB or distilled water is dispensed directly onto the center of column membrane for optimal elution of DNA.

To obtain more concentrated DNA solution, apply 30  $\mu$ l of elution buffer, but the volume lower than 30  $\mu$ l will decrease the yield significantly. Up to 200  $\mu$ l of elution buffer can be applied to mini column, and it will reduce the concentration of DNA.

For larger fragment (>5 kb), use pre-warmed ( $70^{\circ}$ C) elution buffer for best efficiency.

For long-term storage, eluting in Buffer EB (10 mM TrisCl, pH 8.5) or TE, pH 8.0 and storing under  $-20^{\circ}$ C is recommended. When using water for elution, check the pH of water (>pH 7.0) before elution.

#### **B.** Vacuum Protocol

The vacuum pressure should be in the range of the list below. Lower vacuum pressure may reduce DNA yield and purity.

Most commercial vacuum manifold with luer connectors can be used with this protocol.

- I. Attach a Column Type D (mini) to a port of the vacuum manifold tightly.
- 2. Add 10 volumes of Buffer NR to 1 volume of the sample and mix. Transfer the mixture to the mini column by pipetting.

For 50  $\mu$ l reaction, add 500  $\mu$ l of Buffer NR. If the length of DNA is longer than 100 bp, add 5 volumes of Buffer NR.

- 3. Switch on vacuum source to draw the solution through the mini column. When all liquid has been pulled through the mini column, release the vacuum.
- 4. Apply 700  $\mu$ l of Buffer NW and switch on vacuum source. When all liquid has been pulled through the mini column, release the vacuum. Transfer the mini column to a collection tube (provided).
- Centrifuge for an additional I min at full speed to remove residual wash buffer. Transfer the mini column to a new 1.5 ml microcentrifuge tube.

If the mini column has Buffer NW associated with it, centrifuge again for additional I min at full speed before transferring to the new 1.5 ml microcentrifuge tube. Residual ethanol from Buffer NW can inhibit subsequent enzymatic reaction.

## 6. Apply 50 $\mu$ I of Buffer EB or ddH<sub>2</sub>O to the center of the membrane in the mini column, let stand for I min and centrifuge for I min at 10,000 x g above (>12,000 rpm).

Ensure that the Buffer EB or distilled water is dispensed directly onto the center of column membrane for optimal elution of DNA.

To obtain more concentrated DNA solution, apply 30  $\mu$ l of elution buffer, but the volume lower than 30  $\mu$ l will decrease the yield significantly. Up to 200  $\mu$ l of elution buffer can be applied to mini column, and it will reduce the concentration of DNA.

For larger fragment (>5 kb), use pre-warmed (70°C) elution buffer for best efficiency.

For long-term storage, eluting in Buffer EB (10 mM TrisCl, pH 8.5) or TE, pH 8.0 and storing under  $-20^{\circ}$ C is recommended. When using water for elution, check the pH of water (>pH 7.0) before elution.

## **Troubleshooting Guide**

Facts	Possible Causes	Suggestions
Low or no recovery	[Gel SV] Incompletly solubilized gel	The sliced agarose gel should be completely dissolved without any particles. To assist the complete solubilization, mix the tube by vortexing every 2~3 minutes during incubation or increase the incubation time. To use low melt agarose usually results in better recovery.
	[Gel SV] Too high pH of binding mixture	At high pH, the binding of DNA to silica membrane will be significantly reduced. The dye included in Buffer GB indicates the pH of mixture as color change from yellow at optimal pH to brown or purple at abnormally higher pH. If the color of mixture has turned to brown or purple, add $10~\mu l$ of 3 M sodium acetate, (pH 5.0) to the sample and mix. The color of mixture will turn to yellow indicating the correct pH for DNA binding.
	Improper elution buffer	As user's requirement, elution buffer other than Buffer EB can be used. The condition of optimal elution is low salt concentration with alkaline pH (7.0 <ph<9.5). as="" buffer="" conditions.<="" eluent,="" ensure="" or="" other="" th="" that="" used="" was="" water="" when=""></ph<9.5).>
	Elution buffer incorrectly dispensed	Ensure that elution buffer dispensed to the center of membrane. Incorrectly dispensed elution buffer causes inappropriate contact with membrane, followed by poor DNA recovery.
Ligation failure	[Gel SV] Too long or strong exposure to UV on transilluminator	UV destroys the DNA ends. Use UV of long wave length and make the handling time as short as possible when excising the gel slice.

Facts	Possible Causes	Suggestions
Clogged membrane	[Gel SV] Incompletly solubilized gel	See the section 'Incompletely solubilized gel' in the Facts "Low or no recovery"
	[Gel SV] > 1.5% agarose gel is used	For >1.5% agarose gel, 5 volumes of Buffer GB to 1 volume of gel slice should be added. For 100 mg of agarose gel, add 500 $\mu$ l of Buffer GB. If the mini column is clogged, transfer the mixture from the mini column to a 1.5 ml microcentrifuge tube, add 1 volume of Buffer GB to mixture volume. Incubate for 5 minutes at 50°C, proceed again to binding steps.
Enzymatic reaction is not performed well with the	Residual ethanol from Buffer NW remains in eluate	It is essential to remove any residual ethanol included in Buffer NW from column membrane. Centrifuge again for complete removal of ethanol.
purified DNA	Too high salt concentration in eluate	Incubate for 5 minutes after addition of Buffer NW at washing steps.
	Eluate contains denatured ssDNA	For reannealing of ssDNA to dsDNA, incubate ssDNA at 95°C for 2 minutes, and then allow to cool slowly to room temperature.
DNA floats out while loading on agarose gel	Residual ethanol from Buffer NW remains in eluate	It is essential to remove any residual ethanol included in Buffer NW from column membrane. Centrifuge again for complete removal of ethanol.
Non-specific band appears after purification	DNA denatured	Renature the DNA by warming up to 95°C for 1 minute and let cool slowly to room temperature.

#### **Ordering Information**

Products	Scale	Size	Cat. No.	Туре	Products	Scale	Size	Cat. No.	Туре
GeneAll® <b>Hybri</b> o	<b>I-Q<sup>™</sup> f</b> o	r rapid pi	reparation of	plasmid DNA	GeneAll® Exgene	e <sup>TM</sup> for is	olation o	f total DNA	
Plasmid Rapidprep		50	100-150				100	105-101	spin /
	mini	200	100-102	spin		mini	250	105-152	vacuum
					DII CV/		26	105-226	spin /
GeneAll® <i>Expre</i>	<b>b<sup>TM</sup></b> for p	reparatio	n of plasmid l	DNA	Blood SV	Midi	100	105-201	vacuum
		50	101-150	spin /		N443/1	10	105-310	spin /
	mini	200	101-102	vacuum		MAXI	26	105-326	vacuum
		26	101-226				100	106-101	spin /
Plasmid SV	Midi	50	101-250	spin /	C-II () /	mini	250	106-152	vacuum
		100	101-201	vacuum	Cell SV	N400/1	10	106-310	spin /
GeneAll® <i>Exfect</i>	ionTM					MAXI	26	106-326	vacuum
for prepa	aration of	transfect	ion-grade pla	smid DNA			100	108-101	spin /
		50	111-150	spin /		mini	250	108-152	vacuum
Plasmid LE	mini	200	111-102	vacuum	CI: : 0./	- NA: 1:	26	108-226	spin /
(Low Endotoxin)		26	111-226	spin /	Clinic SV	Midi	100	108-201	vacuum
	Midi	100	111-201	vacuum		NANA	10	108-310	spin /
Plasmid EF		20	121-220			MAXI	26	108-326	vacuum
(Endotoxin Free) Midi	100	121-201	spin	Genomic DNA micr	·o	50	118-050	spin	
,							100	117-101	spin /
GeneAll® <i>Expin</i> ¹	M for bur	ification (	of fragment D	NA		mini	250	117-152	vacuum
	' '	50	102-150	spin /	DI . CV /		26	117-226	spin /
Gel SV mini	mini	200	102-102	-F, I Idill 3V	Plant SV	Midi	100	117-201	vacuum
		50	103-150				10	117-310	spin /
PCR SV	mini	200	103-102	vacuum		MAXI	26	117-326	vacuum
		50	113-150	spin /	Soil DNA mini	mini	50	114-150	spin
CleanUp SV	mini	200	113-102	vacuum	Stool DNA mini	mini	50	115-150	spin
		50	112-150		Viral DNA / RNA	mini	50	128-150	spin
					Viral DNA / RNA mini	50	120-130	Shiii	
Combo GP	mini			spin /	FEDE (* DATA		50	138-150	
Combo GP	mini	200	112-102	vacuum	FFPE tissue DNA	mini			spin
Combo GP GeneAll <sup>®</sup> Exgen		200	112-102			mini	50 250	138-150 138-152	spin
	<b>e<sup>TM</sup></b> for is	200	112-102		FFPE tissue DNA  GeneAll® GenEx	mini	50 250	138-150 138-152	spin
		200 olation o	112-102 f total DNA	vacuum	GeneAll® <b>GenE</b> x	mini <sup>TM</sup> for isol	50 250	138-150 138-152	spin thout spin
GeneAll® Exgen	<b>e<sup>TM</sup></b> for is mini	200 olation o	112-102 f total DNA 104-101	vacuum spin /		mini	50 250 lation of	138-150 138-152 total DNA wit	spin thout spin
	<b>e<sup>TM</sup></b> for is	200 olation o 100 250	112-102 f total DNA 104-101 104-152	spin /	GeneAll® <b>GenE</b> x	mini <sup>TM</sup> for isol	50 250 lation of	138-150 138-152 total DNA wit 220-101	spin thout spin solution
GeneAll® Exgen	e <sup>TM</sup> for is mini Midi	200 olation o 100 250 26	112-102 f total DNA 104-101 104-152 104-226	spin / vacuum	GeneAll® <b>GenE</b> x	mini  TM for isol  Sx  Lx	50 250 lation of 1 100 500	138-150 138-152 total DNA wit 220-101 220-105	spin  thout spin  solutior
GeneAll® Exgen	<b>e<sup>TM</sup></b> for is mini	200 olation o 100 250 26 100	112-102 f total DNA 104-101 104-152 104-226 104-201	spin / vacuum spin / vacuum spin / vacuum	GeneAll® <b>GenE</b> x	mini  TM for isol  Sx	50 250 lation of 1 100 500 100	138-150 138-152 total DNA wit 220-101 220-105 220-301	spin thout spin solution
GeneAll® Exgen	e <sup>TM</sup> for is mini Midi MAXI	200 olation o 100 250 26 100 10	112-102 f total DNA 104-101 104-152 104-226 104-201 104-310	spin / vacuum spin / vacuum spin / vacuum	GeneAll® GenEx  GenEx™ Blood	mini  TM for isol  Sx  Lx	50 250 250 250 100 500 100	138-150 138-152 total DNA wit 220-101 220-105 220-301 221-101	spin solutior solutior
GeneAll® Exgen	e <sup>TM</sup> for is mini Midi	200    100   250   26   100   10   26	112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326	spin / vacuum spin / vacuum spin / vacuum spin / vacuum	GenEx <sup>TM</sup> Blood GenEx <sup>TM</sup> Cell	mini  TM for isol  Sx  Lx  Sx  Lx	50 250 250 250 100 500 100 500	138-150 138-152 total DNA wit 220-101 220-105 220-301 221-101 221-105	spin solutior solutior solutior
GeneAll® Exgen Tissue SV	e <sup>TM</sup> for is mini Midi MAXI mini	200    100   250   26   100   10   26   100	112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101	spin / vacuum	GeneAll® GenEx  GenEx™ Blood	mini  TM for isol  Sx  Lx  Sx	50 250 250 100 500 100 100 500	138-150 138-152 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301	spin solutior solutior solutior
GeneAll® Exgen	e <sup>TM</sup> for is mini Midi MAXI	200  olation o  100  250  26  100  10  26  100  250	112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101 109-152	spin / vacuum	GenEx <sup>TM</sup> Blood GenEx <sup>TM</sup> Cell	mini  TM for isol  Sx  Lx  Sx  Lx	50 250 250 100 500 100 100 500 100	138-150 138-152 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301 222-101	spin
GeneAll® Exgen Tissue SV	e <sup>TM</sup> for is mini Midi MAXI mini	200    100   250   26   100   26   100   250   26   26   26   26   26   26   26   2	112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101 109-152 109-226	spin / vacuum	GenEx <sup>TM</sup> Blood GenEx <sup>TM</sup> Cell	mini  TM for isol  Sx  Lx  Sx  Lx  Sx	50 250 250 100 500 100 500 100 100 500	138-150 138-152 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301 222-101 222-105	spin solutior solutior solutior solutior solutior

Products	Scale	Size	Cat. No.	Туре
GeneAll® GenEx	<b>ΓΜ</b> for is	solation of	total DNA	
	Sx	100	227-101	
GenEx <sup>™</sup> Plant	Mx	100	227-201	solution
	Lx	100	227-301	
	Sx	100	228-101	
GenEx <sup>™</sup> Plant plus!	Mx	50	228-250	solution
	Lx	20	228-320	
GeneAll® DirExT/	<sup>M</sup> serie	es		

for preparation of PCR-template without extraction

jo, propordaori o	1 . 0	ipiaco iniciroa	ic onta accion
DirEx <sup>™</sup>	100	250-101	solution
DirEx <sup>™</sup> Fast-Tissue	96 T	260-011	solution
DirEx <sup>™</sup> Fast-Cultured cell	96 T	260-021	solution
DirEx <sup>™</sup> Fast-Whole blood	96 T	260-031	solution
DirEx <sup>™</sup> Fast-Blood stain	96 T	260-041	solution
DirEx <sup>™</sup> Fast-Hair	96 T	260-051	solution
DirEx <sup>™</sup> Fast-Buccal swab	96 T	260-061	solution
DirEx <sup>™</sup> Fast-Cigarette	96 T	260-071	solution

#### GeneAll® RNA series for preperation of total RNA

RiboEx <sup>™</sup>	mini	100	301-001	solution
KIDOEX	mini	200	301-002	SOlution
Hybrid-R <sup>™</sup>	mini	100	305-101	spin
Hybrid-R <sup>™</sup> Blood RNA	4 mini	50	315-150	spin
Hybrid-R <sup>™</sup> miRNA	mini	50	325-150	spin
RiboEx <sup>™</sup> LS	mini	100	302-001	solution
KIDOEX LS	mini	200	302-002	SOIULION
Riboclear™	mini	50	303-150	spin
Riboclear <sup>™</sup> plus!	mini	50	313-150	spin
Ribospin <sup>™</sup>	mini	50	304-150	spin
Ribospin <sup>™</sup> II	mini	50	314-150	
Midospin II		300	314-103	spin
Ribospin <sup>™</sup> vRD	mini	50	302-150	spin
Ribospin ™ vRD <i>plus!</i>	mini	50	312-150	spin
Ribospin <sup>™</sup> vRD II	mini	50	322-150	spin
Ribospin <sup>™</sup> Plant	mini	50	307-150	spin
Ribospin <sup>™</sup> Seed / Fruit	mini	50	317-150	spin
Allspin <sup>™</sup>	mini	50	306-150	spin
RiboSaver™	mini	100	351-001	solution

Products	Scale	Size	C	at. No.	Туре
GeneAll® <b>AmpO</b>	<b>NE<sup>TM</sup></b> for	PCR a	ımţ	olification	
•		250 l	-	501-025	
Taq DNA polymeras	e	500 l	J	501-050	(2.5 U/ <b>µℓ</b> )
		1,000 (	J	501-100	
		250 l	J	502-025	
lpha-Taq DNA polymerase		500 l	J	502-050	(2.5 U/ <b>µl</b> )
		ا 000, ا	J	502-100	
		250 l	J	504-025	
lpha-Pfu DNA polyme	rase	500 l	J	504-050	(2.5 U/ <b>µl</b> )
		1,000 (	J	504-100	
		250 l	J	505-025	
Fast-Pfu DNA polymerase		500 L	J	505-050	(2.5 U/µℓ)
polymerase		1,000 (	IJ	505-100	
		250 l	J	531-025	
Hotstart Taq DNA polymerase		500 L	J	531-050	(2.5 U/ <b>µℓ</b> )
polymerase		1,000 (	J	531-100	
	96 tubes	20 μ	2	521-200	1 120
To a Donation		50 μ(	2	521-500	- lyophilized
Taq Premix		20 μ	2	526-200	1-41
		50 μ <b>(</b>	2	526-500	- solution
	96 tubes	20 μ	2	522-200	h le lii el
or Transfer		50 μ <b>(</b>	2	522-500	- lyophilized
lpha-Taq Premix		20 μ	2	527-200	1
		50 μ(	2	527-500	- solution
HS-Taq Premix	96 tubes	20 μ	2	525-200	1.7
		50 μ(	2	525-500	solution
		20 μ	2	520-200	lyophilized
	96 tubes	50 μ(	2	523-500	solution
Taq Premix (w/o dye)	96 tubes	20 μ	2	524-200	lyophilized
dNTPs mix		500 μ(	2	509-020	2.5 mM eac
dNTPs set (set of dATP, dCTP, dGTP an	d dTTP)	I ml x 4 tube:		509-040	100 mM

Products Scale Size Cat. No. Type	9
-----------------------------------	---

#### **GeneAll<sup>®</sup> AmpMaster<sup>™</sup>** for PCR amplification

Taq Master mix	0.5 ml x 2 tubes	541-010	solution
	0.5 ml x 10 tubes	541-050	solution
lpha-Taq Master mix	0.5 ml x 2 tubes	542-010	solution
	0.5 ml x 10 tubes	542-050	solution
HS-Taq Master mix	0.5 ml x 2 tubes	545-010	solution
	0.5 ml x 10 tubes	545-050	solution
lpha-Pfu Master mix	0.5 ml x 2 tubes	543-010	solution
	0.5 ml x 10 tubes	543-050	solution

#### **GeneAll® HyperScript™** for Reverse Transcription

deficate risperscript for neverse manscription				
Reverse Transcript	ase 10,000 U	601-100	solution	
RT Master mix	$0.5~\mathrm{ml} \times 2~\mathrm{tubes}$	601-710	solution	
RT Master mix with oligo (dT) <sub>20</sub>	$0.5~\mathrm{ml} \times 2~\mathrm{tubes}$	601-730	solution	
RT Master mix with random hexamer	$0.5~\mathrm{ml} \times 2~\mathrm{tubes}$	601-740	solution	
RT Premix	96 tubes, 20 μl	601-602	solution	
RT Premix with oligo (dT) <sub>20</sub>	96 tubes, 20 μ <b>l</b>	601-632	solution	
RT Premix with random hexamer	96 tubes, 20 μl	601-642	solution	
One-step RT-PCR Master mix	$0.5~\mathrm{ml} \times 2~\mathrm{tubes}$	602-110	solution	
One-step RT-PCR Premix	96 tubes, 20 μl	602-102	solution	
First strand Synthesis Kit	50 reaction	605-005	solution	
ZymAll <sup>™</sup> RNase Inhibitor	10,000 ∪	605-010	solution	
ZymAll <sup>TM</sup> RNase Inhibitor	4,000 U	605-004	solution	

#### GeneAll® RealAmp<sup>™</sup> for qPCR amplification

SYBR qPCR Master	200 rxn	20 µl	801-020	solution
mix (2X, Low ROX)	500 rxn	20 µl	801-050	SOIULION
SYBR qPCR Master	200 rxn	20 µl	801-021	1.41
mix (2X, High ROX)	500 rxn	20 µl	801-051	solution

Products Size Cat. No. Type
-----------------------------

#### GeneAll® Protein series

ProtinEx <sup>™</sup> Animal cell / tissue		701-001	solution
PAGESTA <sup>™</sup> Reducing 5X SDS-PAGE Sample Buffer	I mI × I0 tubes	751-001	solution

#### GeneAll $^{\otimes}$ STEAD $\dot{\iota}^{\mathsf{TM}}$ for automatic nucleic acid puritication

	,		,
12 Instrument		GST012	system
24 Instrument		GST024	system
Genomic DNA Cell / Tissue	96	401-104	kit
Genomic DNA Blood	96	402-105	kit
Total RNA	96	404-304	kit
Viral DNA / RNA	96	405-322	kit
CFC Seed DNA / RNA	96	406-C02	kit
Genomic DNA Plant	96	407-117	kit
Soil DNA	96	408-114	kit

## NOTE



#### www.geneall.com

GeneAll Bldg., 303-7 Dongnam-ro, Songpa-gu, Seoul, South Korea 138-859 E-mail: sales@geneall.com

Tel: 82-2-407-0096 Fax: 82-2-407-0779

©2017 GeneAll Biotechnology, All right reserved