

Ver 2.0

HBI150

g

Handbook for

■ Soil DNA mini

exgene™

DNA PURIFICATION HANDBOOK


GeneAll

Customer & Technical Support

Do not hesitate to ask us any question.

We thank you for any comment or advice.

Contact us at

www.geneall.com

Tel : 82-2-407-0096

Fax : 82-2-407-0779

E-mail(Order/Sales) : sales@geneall.com

E-mail(Tech. Info.) : tech@geneall.com

Visit GeneAll® Community

www.geneall.com

www.geneall.co.kr

Trademarks

Hybrid-Q™, Exprep™, Exfection™, Expin™, Exgene™, GenEx™, DirEx™, Hybrid-R™, RiboEx™, Ribospin™, Riboclear™, Allspin™, RiboSaver™, EzClear™, EzSep™, EzPure™, EzPass™, AmpONE™, AmpMaster™, RealAmp™, HyperScript™, ProtinEx™, PAGESTA™, STEADi™ are trademarks of GeneAll Biotechnology co., Ltd.

© 2016 GeneAll Biotechnology, all right reserved.

This protocol handbook is included in :

GeneAll® Exgene™ Soil DNA mini (I14-I50)

Visit www.geneall.com or www.geneall.co.kr for FAQ, QnA and more information.

Sample pulverization step

Add up to 500 mg of soil sample to a Powerbead™ tube.
Add 550 ul of buffer SL.
Pulverize the sample.
Centrifuge at $\geq 10,000 \times g$ for 10 minutes.

Inhibitor removal step

Transfer the supernatant to a 1.5 ml microcentrifuge tube.
Add 50 ul of buffer RH.
Add 300 ul of buffer PD and mix well.
Centrifuge at $\geq 10,000 \times g$ for 5 minutes.

DNA binding step

Transfer the supernatant to a 2.0 ml microcentrifuge tube.
Add 900 ul of buffer TB.
Apply the mixture into a mini spin column and centrifuge at $\geq 10,000 \times g$ for 30 seconds.


Washing step

Add 500 ul of buffer NW and
Centrifuge at $\geq 10,000 \times g$ for 30 seconds.
Centrifuge at $\geq 10,000 \times g$ for 1 minute.

DNA elution

Add ~50 ul of buffer EB to the center of the membrane.
Centrifuge at $\geq 10,000 \times g$ for 1 minute.

Brief Protocol



INDEX

Index	05
Kit Contents	06
Materials Not Provided	
Quality Control	07
Storage Conditions	
Precautions	
Product Disclaimer	
Product Specifications	08
Product Description	09
Protocol for Exgene™ Soil DNA mini	10
Troubleshooting Guide	12

KIT CONTENTS

Components	Quantity	Storage
Buffer SL	30 ml	Room temperature
Buffer RH	3 ml	
Buffer PD	17 ml	
Buffer TB	50 ml	
Buffer NW (concentrate) * †	6 ml	
Buffer EB	15 ml	
Powerbead™ tube	50	
GeneAll® Column type G (with collection tube)	50	
1.5 ml microcentrifuge tube	100	
2.0 ml microcentrifuge tube	50	

* Before using for the first time, add absolute ethanol (ACS grade or better) into buffer NW as indicated on the bottle.

† Contains sodium azide as a preservative.

MATERIALS NOT PROVIDED

Reagent

- Absolute ethanol, ACS grade or better

Disposable material

- Pipet tips
- Disposable gloves

Equipment

- Precellys® 24 (Bertin, France) equipment or any equivalent
- Microcentrifuge
- Suitable protector (ex; lab coat, disposable gloves, goggles, etc)

QUALITY CONTROL

GeneAll® Exgene™ Soil DNA mini is manufactured in strictly clean condition, and its degree of cleanness is monitored periodically. For consistency of product, the quality certification process is carried out from lot to lot thoroughly and only the qualified is approved to be delivered.

STORAGE CONDITIONS

GeneAll® Exgene™ Soil DNA mini should be stored at room temperature (15~25°C). But prolonged storage at high temperature over 30°C can reduce the performance of the kit.

In cold ambient condition, buffer RH and TB may exhibit salt precipitation and this will cause reduction of DNA recover-yields. If so, heat the bottle with occasional swirling in 37°C water bath until completely dissolved.

All components are stable for 1 year.

Keep out of direct sunlight.

PRECAUTIONS

The buffers included in GeneAll® Exgene™ Soil DNA mini contain irritant which is harmful when in contact with skin or eyes, or when inhaled or swallowed. Care should be taken during handling. Always wear gloves and eye protector, and follow standard safety precautions. In case of contact, wash immediately with plenty of water and seek medical advice.

Buffer TB contains chaotropes. It can form highly reactive compounds when combined with bleach. Do NOT add bleach or acidic solutions directly to the sample-preparation waste.

PRODUCT DISCLAIMER

GeneAll® Exgene™ Soil DNA mini is for research use only, not for use in diagnostic procedure.

Product Specifications

Specification	Exgene™ Soil DNA mini
Type	Spin
Maximum amount of starting samples	500 mg soil sample
Maximum loading volume of spin column	700 ul
Minimum elution volume	30 ul
Maximum binding capacity	100 ug

Product Description

GeneAll® Exgene™ Soil DNA mini provides a convenient method for the isolation of total DNA from soil samples. This kit utilizes the powerful beads, the optimized buffer system and the advanced silica binding technology to purify nucleic acid suitable for many applications. These complex systems of this kit can deal with a number of different types of samples in the soil including plant tissues, bacteria, fungi spores and others. Also, it removes a humic acid and other PCR inhibitors from various soil samples efficiently. The humic acid, which is a sort of brownish colour, is a critical factor for soil treating experiments and if remained in eluate, this can have a negative effect on the DNA downstream applications.

GeneAll® Exgene™ Soil DNA mini provides a tube including powerful beads for strong pulverization. Soil samples are placed in this tube with lysis buffer, buffer SL, and crushed by bead-beater or vortex. After centrifugation, supernatant is mixed with precipitation buffer, buffer RH and buffer PD, to precipitate humic acid and protein. Then, the separated DNA part, supernatant, blend into the binding buffer, buffer TB, and DNA is bound on the silica membrane through centrifugation. Following washing step with buffer NW, the bound DNA is eluted by buffer EB. Purified DNA can be directly applicable in conventional PCR, restriction analysis, electrophoresis, and any other downstream applications.

PROTOCOL FOR

Exgene™ Soil DNA mini

- 1. Add up to 500 mg of soil sample to a Powerbead™ tube.**
- 2. Add 550 ul of buffer SL to the tube.**
- 3. Homogenize the sample in the Precellys® 24 (Bertin, France) equipment for twice of 23 seconds at 6500 rpm.**

Alternatively, secure tubes horizontally on a flat-bed vortex pad with tape and vortex at maximum speed for 10 minutes.
- 4. Centrifuge at $\geq 10,000 \times g$ for 10 minutes at room temperature and carefully transfer the supernatant to a 1.5 ml microcentrifuge tube (provided).**
- 5. Add 50 ul of buffer RH.**
- 6. Add 300 ul of buffer PD and mix well by vortexing.**
- 7. Centrifuge at $\geq 10,000 \times g$ for 5 minutes at room temperature and carefully transfer the supernatant to a 2.0 ml microcentrifuge tube (provided).**

Small pellet containing humic acid, cell debris, and protein can be formed in the collection tube after centrifugation. Be careful not to disturb this pellet.
- 8. Add 900 ul of buffer TB and mix well by vortexing.**

If buffer TB precipitation, pre-heat in a 56°C water bath to dissolve completely.
- 9. Transfer up to 700 ul of the mixture to a mini spin column.**
- 10. Centrifuge at $\geq 10,000 \times g$ for 30 seconds at room temperature.**

Discard the pass-through and reinsert the mini spin column back into the same tube.

11. Repeat two more times step 9~10 using the remainder of the sample.

12. Add 500 ul of buffer NW to the mini spin column.

13. Centrifuge at $\geq 10,000 \times g$ for 30 seconds at room temperature.

Discard the pass-through and reinsert the mini spin column back into the same tube.

14. Centrifuge at maximum speed for 1 minute at room temperature to remove residual wash buffer.

Transfer the mini spin column to a new 1.5 ml microcentrifuge tube (provided).

Residual ethanol may interfere with downstream reactions. Care must be taken at this step for eliminating the carryover of buffer NW.

15. Add 50 ul of buffer EB to the center of the membrane in the mini spin column.

Incubate for 1 minute at room temperature. Centrifuge at $\geq 10,000 \times g$ for 1 minute at room temperature.

Elution volume can be decreased to 30 ul for high concentration of DNA, but this will slightly decrease in overall DNA yield. If maximum recovery of DNA is preferred or the starting materials contain large amount of DNA, elution can be done in 200 ul of buffer EB.

Troubleshooting Guide

Facts	Possible Causes	Suggestions
Low or no recovery	Too much starting material	Too much starting material lead to inefficient lysis, followed by poor DNA yields. Reduce the amount of starting material.
	Insufficient Homogenization	Check the step 3 of protocol. Insufficient homogenization time and condition is related to low recovery yield.
Low efficiency of DNA amplification	Excess amount of template DNA	An excess amount of template DNA will inhibit a PCR reaction. The template DNA is needed to dilute.
Eluate does not perform well in the downstream application	Residual ethanol remains in eluate	To remove any residual ethanol included in buffer NW from mini spin column membrane, centrifuge again for complete removal of ethanol.
DNA eluate is brown	Humic acid is not be removed completely	With certain samples, a little humic acid can be remained in the eluate. In this case, we recommend using a GeneAII® Expin™ CleanUp SV kit to purify contaminated eluate.

Ordering Information

Products	Scale	Size	Cat. No.	Type	Products	Scale	Size	Cat. No.	Type	
GeneAll® Hybrid-Q™ for rapid preparation of plasmid DNA					GeneAll® Exgene™ for isolation of total DNA					
Plasmid Rapidprep	mini	50	100-150	spin	Blood SV	mini	100	105-101	spin / vacuum	
		200	100-102				250	105-152		
GeneAll® Exprep™ for preparation of plasmid DNA	mini	50	101-150	spin / vacuum		Midi	26	105-226	spin / vacuum	
		200	101-102				100	105-201		
	Plasmid SV	Midi	26	101-226		spin / vacuum	MAXI	10	105-310	spin / vacuum
			50	101-250				26	105-326	
		100	101-201							
GeneAll® Exfection™ for preparation of transfection-grade plasmid DNA					GeneAll® GenEx™ for isolation of total DNA without spin column					
Plasmid LE (Low Endotoxin)	mini	50	111-150	spin / vacuum	Cell SV	mini	100	106-101	spin / vacuum	
		200	111-102				250	106-152		
	Midi	26	111-226	spin / vacuum		MAXI	10	106-310	spin / vacuum	
100		111-201	26				106-326			
Plasmid EF (Endotoxin Free)	Midi	20	121-220	spin		mini	100	108-101	spin / vacuum	
		100	121-201				250	108-152		
GeneAll® Expin™ for purification of fragment DNA	mini	50	111-150	spin / vacuum	Midi	Midi	26	108-226	spin / vacuum	
		200	102-150				100	108-201		
	Gel SV	mini	50	103-150	spin / vacuum	MAXI	10	108-310	spin / vacuum	
200			103-102	26			108-326			
PCR SV	mini	50	113-150	spin / vacuum	Clinic SV	mini	100	108-101	spin / vacuum	
		200	113-102				250	108-152		
CleanUp SV	mini	50	112-150	spin / vacuum		Midi	26	108-226	spin / vacuum	
		200	112-102				100	108-201		
Combo GP	mini	50	117-150	spin / vacuum		MAXI	10	108-310	spin / vacuum	
		200	117-102				26	108-326		
GeneAll® Exgene™ for isolation of total DNA					GeneAll® GenEx™ for isolation of total DNA without spin column					
Tissue SV	mini	100	104-101	spin / vacuum	Plant SV	mini	100	114-150	spin	
		250	104-152				250	117-152		vacuum
	Midi	26	104-226	spin / vacuum		Midi	26	117-226	spin / vacuum	
		100	104-201				100	117-201		vacuum
MAXI	10	104-310	spin / vacuum	MAXI		10	117-310	spin / vacuum		
	26	104-326				26	117-326		vacuum	
Tissue plus! SV	mini	100	109-101	spin / vacuum	Soil DNA mini	mini	50	115-150	spin	
		250	109-152				50	128-150		spin
	Midi	26	109-226	spin / vacuum	Stool DNA mini	mini	50	138-150	spin	
		100	109-201				50	138-150		spin
	MAXI	10	109-310	spin / vacuum	Viral DNA / RNA	mini	50	138-150	spin	
		26	109-326				250	138-152		spin

Products	Scale	Size	Cat. No.	Type
GeneAll® GenEx™ for isolation of total DNA				
GenEx™ Plant	Sx	100	227-101	solution
	Mx	100	227-201	
	Lx	100	227-301	
GenEx™ Plant plus!	Sx	100	228-101	solution
	Mx	50	228-250	
	Lx	20	228-320	

GeneAll® DirEx™ series
for preparation of PCR-template without extraction

DirEx™		100	250-101	solution
DirEx™ Fast-Tissue		96 T	260-011	solution
DirEx™ Fast-Cultured cell		96 T	260-021	solution
DirEx™ Fast-Whole blood		96 T	260-031	solution
DirEx™ Fast-Blood stain		96 T	260-041	solution
DirEx™ Fast-Hair		96 T	260-051	solution
DirEx™ Fast-Buccal swab		96 T	260-061	solution
DirEx™ Fast-Cigarette		96 T	260-071	solution

GeneAll® RNA series for preparation of total RNA

RiboEx™	mini	100	301-001	solution
		200	301-002	
Hybrid-R™	mini	100	305-101	spin
Hybrid-R™ Blood RNA mini		50	315-150	spin
Hybrid-R™ miRNA	mini	50	325-150	spin
RiboEx™ LS	mini	100	302-001	solution
		200	302-002	
Riboclear™	mini	50	303-150	spin
Riboclear™ plus!	mini	50	313-150	spin
Ribospin™	mini	50	304-150	spin
Ribospin™ II	mini	50	314-150	spin
		300	314-103	
Ribospin™ vRD	mini	50	302-150	spin
Ribospin™ vRD plus!	mini	50	312-150	spin
Ribospin™ vRD II	mini	50	322-150	spin
Ribospin™ Plant	mini	50	307-150	spin
Ribospin™ Seed / Fruit	mini	50	317-150	spin
Allspin™	mini	50	306-150	spin
RiboSaver™	mini	100	351-001	solution

Products	Scale	Size	Cat. No.	Type
GeneAll® AmpONE™ for PCR amplification				
Taq DNA polymerase		250 U	501-025	(2.5 U/μl)
		500 U	501-050	
		1,000 U	501-100	
α-Taq DNA polymerase		250 U	502-025	(2.5 U/μl)
		500 U	502-050	
		1,000 U	502-100	
α-Pfu DNA polymerase		250 U	504-025	(2.5 U/μl)
		500 U	504-050	
		1,000 U	504-100	
Fast-Pfu DNA polymerase		250 U	505-025	(2.5 U/μl)
		500 U	505-050	
		1,000 U	505-100	
Hotstart Taq DNA polymerase		250 U	531-025	(2.5 U/μl)
		500 U	531-050	
		1,000 U	531-100	
Taq Premix	96 tubes	20 μl	521-200	lyophilized
		50 μl	521-500	
		20 μl	526-200	
50 μl	526-500			
α-Taq Premix	96 tubes	20 μl	522-200	lyophilized
		50 μl	522-500	
		20 μl	527-200	
50 μl	527-500			
HS-Taq Premix	96 tubes	20 μl	525-200	solution
		50 μl	525-500	
α-Pfu Premix	96 tubes	50 μl	523-500	solution
		20 μl	524-200	
Taq Premix (w/o dye)	96 tubes	20 μl	524-200	lyophilized
dNTPs mix		500 μl	509-020	2.5 mM each
dNTPs set (set of dATP, dCTP, dGTP and dTTP)		1 ml x 4 tubes	509-040	100 mM

Products	Scale	Size	Cat. No.	Type
----------	-------	------	----------	------

GeneAll® AmpMaster™ for PCR amplification

Taq Master mix	0.5 ml x 2 tubes	541-010	solution
	0.5 ml x 10 tubes	541-050	solution
α-Taq Master mix	0.5 ml x 2 tubes	542-010	solution
	0.5 ml x 10 tubes	542-050	solution
HS-Taq Master mix	0.5 ml x 2 tubes	545-010	solution
	0.5 ml x 10 tubes	545-050	solution
α-Pfu Master mix	0.5 ml x 2 tubes	543-010	solution
	0.5 ml x 10 tubes	543-050	solution

GeneAll® HyperScript™ for Reverse Transcription

Reverse Transcriptase	10,000 U	601-100	solution
RT Master mix	0.5 ml x 2 tubes	601-710	solution
RT Master mix with oligo (dT) ₂₀	0.5 ml x 2 tubes	601-730	solution
RT Master mix with random hexamer	0.5 ml x 2 tubes	601-740	solution
RT Premix	96 tubes, 20 μl	601-602	solution
RT Premix with oligo (dT) ₂₀	96 tubes, 20 μl	601-632	solution
RT Premix with random hexamer	96 tubes, 20 μl	601-642	solution
One-step RT-PCR Master mix	0.5 ml x 2 tubes	602-110	solution
One-step RT-PCR Premix	96 tubes, 20 μl	602-102	solution
First strand Synthesis Kit	50 reaction	605-005	solution
ZymAll™ RNase Inhibitor	10,000 U	605-010	solution
ZymAll™ RNase Inhibitor	4,000 U	605-004	solution

GeneAll® RealAmp™ for qPCR amplification

SYBR qPCR Master mix (2X, Low ROX)	200 rxn	20 μl	801-020	solution
	500 rxn	20 μl	801-050	
SYBR qPCR Master mix (2X, High ROX)	200 rxn	20 μl	801-021	solution
	500 rxn	20 μl	801-051	

Products	Size	Cat. No.
----------	------	----------

GeneAll® Protein series

ProteinEx™ Animal cell / tissue	100 ml	701-001	solution
PAGESTA™ Reducing 5X SDS-PAGE Sample Buffer	1 ml x 10 tubes	751-001	solution

GeneAll® STEADi™ for automatic nucleic acid purification

STEADi™ 12 Instrument		GST012
STEADi™ 24 Instrument		GST024
STEADi™ Genomic DNA Cell / Tissue Kit	96	401-104
STEADi™ Genomic DNA Blood Kit	96	402-105
STEADi™ Bacteria DNA Kit	96	403-106
STEADi™ Total RNA Kit	96	404-304
STEADi™ Viral DNA / RNA Kit	96	405-322
STEADi™ CFC Seed DNA / RNA Kit	96	406-C02



GeneAll

GENEALL BIOTECHNOLOGY CO., LTD

www.geneall.com

GeneAll Bldg., 303-7 Dongnam-ro,
Songpa-gu, Seoul, South Korea | 38-859
E-mail : sales@geneall.com

Tel : 82-2-407-0096
Fax : 82-2-407-0779

©2016 GeneAll Biotechnology, All right reserved

Edited by **BnP**
Designed by **KDY**
2016.08